


B Product Specifications



Model#	NL-100
Type	Breast Marker
Marking Edges	38mm & 42mm
Materials	Polycarbonate

C Sterilization and Disinfection

Sterilization and High-Level Disinfection Recommendations:

Pasteurization: Pasteurize at 160°F-170°F (71°C-77°C) for no less than 30 minutes.

Gamma Radiation: Irradiate according to validated parameters. Some discoloration may occur.

Ethylene Oxide: Sterilize according to validated parameters. Do not exceed 131°F (55°C). Allow ample aeration time in a well-ventilated area to dissipate the absorbed gas.

Chemical Disinfectants: Disinfect according to validated parameters. Follow the chemical manufacturer's recommendation for temperature and soak time. Chemical disinfection should be followed by a sterile water rinse.

Gas Plasma: Sterilize according to validated parameters.



Instrumentation

Industries, Inc.

2990 Industrial Boulevard • Bethel Park, PA 15102

U.S. Toll Free: 1-800-633-8577

Business: 412-854-1133 • **Fax:** 412-854-5668

E-mail: iiimed@iiirespiratory.com

www.iiirespiratory.com

ECO 1209 Rev. A



Instrumentation
Industries, Inc.

Directions for Use

Breast Marker

NL-100



A

NL-100 Breast Marker Directions for Use

Fig. 1





Fig. 2



Fig. 3



NL-100 Breast Markers are:

- Clean, but **NOT** sterile.
- Single-Patient Use **ONLY**.
- Packaged in quantities of 25.

Directions for use of the **NL-100 Breast Marker**:

1. **STERILIZE THE MARKER USING THE RECOMMENDATIONS IN SECTION C.**
2. The **NL-100 Breast Marker** provides two marking edges; 42mm and 38mm
3. Choose the marking size with or without methylene blue as a marking agent.
4. Place the **NL-100 Breast Marker** over the nipple areola complex, centering through the cross hairs. *Fig. 1*
5. Apply firm downward pressure and rotate the **NL-100 Breast Marker** in a clockwise and counter clockwise manner. *Fig. 2*
6. Remove the **NL-100 Breast Marker**. The areola stencil can now be visualized. *Fig. 3*