


A RTC 24-V / RTC 24-V KIT — Installation Directions

Note: All illustrations are examples.
Consult facility validated procedures for desired placement.

RTC 24-V — Placed at Endotracheal Tube


1. Connect RTC 24-V Inlet to Patient Wye, as shown. ①
2. Connect RTC 24-V Outlet to Endotracheal Tube. ②
3. Ensure that Rubber Cap is placed on MDI Activator Port until MDI canister is inserted. ③



Direction arrow on RTC 24-V must be pointing in the direction of airflow.

RTC 24-V KIT — Placed at Endotracheal Tube — RTC 24-V MDI Adapter • KC 004-5 Tubing • KC 2124 Adapter


1. Assemble RTC 24-V, KC 004-5, and KC 2124 before connecting to the Patient Circuit.
 - a. Connect RTC 24-V Outlet to KC 004-5. ①
 - b. Connect RTC 24-V/KC 004-5 to KC 2124. ②
2. Connect RTC 24-V/KC 004-5/KC 2124 Assembly to Patient Wye at RTC 24-V Inlet, as shown. ③
3. Connect RTC 24-V/KC 004-5/KC 2124 Assembly to Endotracheal Tube at KC 2124. ④
4. Ensure that Rubber Cap is placed on MDI Activator Port until MDI canister is inserted. ⑤



Direction arrow on RTC 24-V must be pointing in the direction of airflow.

RTC 24-V KIT — Placed Inline — RTC 24-V MDI Adapter • KC 004-5 Tubing • KC 2124 Adapter

1. Assemble RTC 24-V, KC 004-5, and KC 2124 before connecting to the Patient Circuit.
 - a. Connect KC 2124 to RTC 24-V Inlet. ①
 - b. Connect KC 2124/RTC 24-V to KC 004-5 at RTC 24-V Outlet. ②
2. Connect KC 2124/RTC 24-V/KC 004-5 Assembly to Inspiratory Limb of Patient Circuit at KC 2124. ③
3. Connect KC 2124/RTC 24-V/KC 004-5 Assembly to Patient Wye at KC 004-5. ④
4. Ensure that Rubber Cap is placed on MDI Activator Port until MDI canister is inserted. ⑤




Direction arrow on RTC 24-V must be pointing in the direction of airflow.

B Canister Insertion

Standard Canister

1. Install RTC 24-V in Patient Circuit. (panel A)
2. Remove Rubber Cap from RTC 24-V.
3. Insert Pressurized Medicine Canister into top of RTC 24-V.
4. Ensure that Canister Activator is inserted into Canister Activator Port.




Note: Standard Canisters are not stable after insertion into RTC 24-V.

Constant pressure **must** be applied between Standard Canister and RTC 24-V.

Canister with Dose Counter

1. Install RTC 24-V in Patient Circuit. (panel A)
2. Remove Rubber Cap from RTC 24-V.
3. Insert Pressurized Medicine Canister into top of RTC 24-V with Dose Counter facing away from Direction Arrow.


4. Ensure that Canister Activator is inserted into Canister Activator Port.
5. Ensure that Dose Counter Activator is inserted into Dose Counter Activator Slot.



C Canister Activation

- Prior to activation, make sure that the Pressurized Medicine Canister is fully seated in RTC 24-V.
- While activating the Pressurized Medicine Canister, press the canister straight down without tilting.
- Two hands may be necessary to keep canister aligned properly.

1. At the beginning of the inspiratory cycle, depress the Pressurized Medicine Canister by placing fingers underneath and pressing top of canister with thumb. If canister has a Dose Counter, it will be tripped by the Dose Counter Activator.



2. Remove Canister from RTC 24-V.
3. Replace Rubber Cap.

D Indications for Use

RTC 24-V — Metered Dose Inhaler Adapter

RTC 24-V KIT — Metered Dose Inhaler Adapter Kit

The Instrumentation Industries, Inc. RTC 24-V MDI Adapter and RTC 24-V KIT MDI Adapter Kit are actuators for intermittent delivery of prescribed aerosol medication dispensed in cylindrical-style metered dose inhalers. The RTC 24-V MDI Adapter and RTC 24-V KIT MDI Adapter Kit are intended for use only when connected to ventilator tubing or tracheal tubes.

The RTC 24-V MDI Adapter and RTC 24-V KIT MDI Adapter Kit are intended to be prescribed for any patient who is ventilator-dependent and to whom a metered dose inhaler has been prescribed. The expected clinical environment for the RTC 24-V MDI Adapter and RTC 24-V KIT MDI Adapter Kit is Critical Care and/or long term or short term ventilation.

The components of the RTC 24-V KIT consist of one RTC 24-V MDI Adapter, one four inch length of 22 mm inner-diameter tubing, and one 15 mm inner-diameter/22 mm outer-diameter adapter.

The RTC 24-V MDI Adapter and RTC 24-V KIT MDI Adapter Kit are intended for single patient reuse.

These devices are intended for sale by or on the order of a physician.

E Sterilization/Disinfection

These products are designed for **Single Patient Reuse**. (Can be used to repeatedly administer medications to one single patient while installed into a breathing circuit)

They are packaged as **non-sterile**; Instrumentation Industries, Inc. recommends processing before use. After installation into a breathing circuit, cleaning is not necessary. Discard when replacing the circuit according to facility guidelines.

Recommendations for RTC 24-V and RTC 24-V KIT:

Carefully follow facility guidelines for preparing products and solutions for the sterilizing and disinfecting methods and equipment being used.

Sterilization and High Level Disinfection Procedures have been validated for the following equipment and methods:

Sterilization

STERRAD 100S

Full short sterilization cycle (Hydrogen peroxide gas plasma, 55 minutes, 45°C - 55°C).

High Level Disinfection

Pasteurization

Pasteurize at 70°C +/- 3°C (153° - 163°F) for a minimum of 30 minutes.


Cidex® OPA

Prepare Cidex® OPA solution. Use a syringe to fill all crevices with Cidex® OPA, then totally submerge adapters in the disinfectant. Follow the manufacturer's instructions for temperature and soak time requirements.

Always carefully examine the physical integrity of this product after processing. Do not use if cracking or crazing is seen, or if parts fit together improperly after processing.

F Product Specifications



RTC 24-V



Model #	RTC 24-V
Type	Metered Dose Inhaler Adapter
Use	Single Patient Reuse
Port Sizes	
Inlet	22mm I.D.
Outlet	15mm I.D./22mm O.D.
Materials	
Body	Polystyrene-Butadiene
Cap	Thermoplastic Rubber
Color	
Body	Clear
Cap	Clear

RTC 24-V KIT Contains one of each

Model #	RTC 24-V	See Above
Model #	KC 004-5	
Type	4" Tubing	
Use	Single Patient Reuse	
Size	22mm I.D.	
Material	EVA Copolymer	
Color	Transparent	
Model #	KC 2124	
Type	Adult Adapter	
Use	Single Patient Reuse	
Port Sizes	15mm I.D./22mm O.D. x 15mm I.D./22mm O.D.	
Material	Polystyrene-Butadiene	
Color	Clear	



 Instrumentation Industries, Inc.

MDI Adapter

Installation & Usage Directions



RTC 24-V RTC 24-V KIT

For use with metered dose inhaler canisters, with or without a dose counter, with standard size metal (.109 inch diameter) or plastic (.124 inch diameter) stems.

Caution: Federal law restricts this device to sale by or on the order of a physician.

 Instrumentation Industries, Inc.

2990 Industrial Blvd.
Bethel Park, PA 15102
US Toll Free: 1-800-633-8577
Business: 1-412-854-1133
Fax: 1-877-633-8661

www.iiimedical.com