


F Cleaning Recommendations

Note: Do not clean Reservoir Bag with alcohol.
Reservoir bag is Single Patient Use.

1. Disassemble the SG 066 O₂ Enrichment Kit as illustrated.



2. Wash the SG 066-1 O₂ Manifold, the KC 024 Tubing, and the SG 193-M 22mm O.D. x 30mm I.D. Adapter in warm soapy water. If necessary, use soft, long-bristle bottle brushes to remove all foreign matter.
Caution: Do not attempt to clean the interior of the Reservoir Bag; if using brushes, be sure not to damage the Ambient Air Inlet Port on the O₂ Manifold.
3. Rinse the cleaned parts with water, making certain all soap residue is removed.
4. Disinfect the SG 066-1 O₂ Manifold, the KC 024 Tubing, and the SG 193-M 22mm O.D. x 30mm I.D. Adapter using one of the following methods:
 - a. Soak for 30 minutes in a solution of one part white vinegar to three parts water.
Rinse thoroughly.
 - b. Soak in an approved chemical disinfectant, following the manufacturer's instructions.
Rinse thoroughly.
5. Shake excess moisture from all parts and place on a clean towel to air dry. Do not wipe dry with a towel; do not use a blow dryer. All parts must be completely dry before attaching to the ventilator.

G Product Specifications

	Model #	SG 066
	Type	Reusable O ₂ Enrichment Kit (Single Patient Use Bag)
	Bag Volume	1-Liter
	Materials	Polysulfone Silicone Rubber Anodized Aluminum Neoprene Bag
Replacement Parts		
	Model #	SG 066-2N
	Type	Single Patient Use Breathing Bag
	Volume	1-Liter
	Model #	KC 024
	Type	Reusable KC Tubing
	Length	24"
	Model #	SG 066-1
	Type	Reusable O ₂ Manifold
	Materials	Polysulfone Silicone Rubber
	Model #	SG 193-M
	Type	Reusable Adapter
	Size	22mm O.D. x 30mm I.D.
	Materials	Polysulfone



Visit iiimedical.com/symbols.pdf for the Glossary of Symbols used in Instrumentation Industries, Inc. labeling.

Instrumentation Industries, Inc.

2990 Industrial Blvd. • Bethel Park, PA 15102

US Toll Free: 1-800-633-8577

Business: 1-412-854-1133

US Toll Free Fax: 1-877-633-8661

Fax: 1-412-854-5668

E-mail: sales@iiimedical.com

www.iiimedical.com

 Instrumentation Industries, Inc.

O₂ Enrichment Kit

Installation & Usage Directions

SG 066



Reusable

— Not made with Natural Rubber Latex

— Not made with Di(2-ethylhexyl) phthalate (DEHP)

Made in USA!

1-23

Rx ONLY

NON STERILE

MR UNSAFE

ECO 2783 Rev. E

A Indications for Use

The SG 066 Oxygen (O₂) Enrichment Kit is designed for use on certain ventilators to provide higher concentrations of oxygen to the patient by mixing room air with an oxygen source at the proximal airway.

B Contraindications

None known.

C Cautions & Warnings

Warning: Never operate the SG 066 O₂ Enrichment Kit without securing it to a ventilator.


Warning: The SG 066 O₂ Enrichment Kit is designed to operate with a hospital grade oxygen supply. Pressure should not exceed 10 psi at the O₂ Manifold Inlet. Oxygen flow to the SG 066 O₂ Enrichment Kit should not exceed 10 liters per minute.

Caution: Ensure that the Intake Port is correctly located, so that it cannot be blocked.

D SG 066 O₂ Enrichment Kit Directions for Set-Up

Directions for Set-up:

- Slide the SG 066 O₂ Enrichment Kit onto the mounting rail using the SG 066-3 O₂ Manifold Bracket. ① Ensure that the Ambient Air Intake Port is facing toward the front of the ventilator.
- Connect the SG 193-M Adapter to the KC 024 Tubing Outlet. ②
- Connect Small Bore Tubing (not included) to a medical grade, pressure compensated O₂ Flowmeter (not included). ③
- Connect the other end of the Small Bore Tubing to the O₂ Inlet at the top of the O₂ Manifold. ④
- Connect the KC 024 Tubing/SG 193-M Adapter Assembly to Inlet Filter on Ventilator. ⑤



E SG 066 O₂ Enrichment Kit Directions for Operation

A physician must determine if there is a need for supplemental oxygen.

The SG 066 O₂ Enrichment Kit can deliver up to 100% oxygen concentration at the proximal airway.



- Determine the Minute Volume using the following formula:

$$\dot{V}_e = (\text{BPM})(V_t)$$

\dot{V}_e = Minute Volume in liters/min-
 V_t = Tidal Volume
 BPM = Breaths Per Minute

- Determine the Inspiratory:Expiratory (I:E) ratio (1:1 or 1:2).
- After determining the I:E ratio, use the tables below:
 - Locate the percentage of oxygen required on the Y-axis.
 - Locate the predetermined minute volume on the diagram.
 - Determine the flow rate of oxygen needed to achieve the desired oxygen concentration by locating the point on the X-axis directly below the established intersection.
- Analyze the delivered oxygen concentration for accuracy.

Note: These tables are not applicable if an oxygen concentrator is in use.



Note: During normal operation of the SG 066 O₂ Enrichment Kit, the reservoir bag will partially collapse with each cycle. This is due to the method at which the oxygen is blended with the room air.