


## F Cleaning Recommendations

**Note: Do not clean Reservoir Bag with alcohol. Reservoir bag is Single Patient Use.**

1. Disassemble the SG 066 O<sub>2</sub> Enrichment Kit as illustrated.



2. Wash the SG 066-1 O<sub>2</sub> Manifold, the KC 024 Tubing, and the SG 193-M 22mm O.D. x 30mm I.D. Adapter in warm soapy water. If necessary, use soft, long-bristle bottle brushes to remove all foreign matter.  
**Caution: Do not attempt to clean the interior of the Reservoir Bag; if using brushes, be sure not to damage the Ambient Air Inlet Port on the O<sub>2</sub> Manifold.**
3. Rinse the cleaned parts with water, making certain all soap residue is removed.
4. Disinfect the SG 066-1 O<sub>2</sub> Manifold, the KC 024 Tubing, and the SG 193-M 22mm O.D. x 30mm I.D. Adapter using one of the following methods:
  - a. Soak for 30 minutes in a solution of one part white vinegar to three parts water.  
**Rinse thoroughly.**
  - b. Soak in an approved chemical disinfectant, following the manufacturer's instructions.  
**Rinse thoroughly.**
5. Shake excess moisture from all parts and place on a clean towel to air dry. **Do not wipe dry with a towel; do not use a blow dryer.** All parts must be completely dry before attaching to the ventilator.

## G Product Specifications

	Model #	SG 066
	Type	Reusable O <sub>2</sub> Enrichment Kit (Single Patient Use Bag)
	Bag Volume	1-Liter
	Materials	Polysulfone Silicone Rubber Anodized Aluminum Neoprene Bag
Replacement Parts		
	Model #	SG 066-2N
	Type	Single Patient Use Breathing Bag
	Volume	1-Liter
	Model #	KC 024
	Type	Reusable KC Tubing
	Length	24"
	Model #	SG 066-1
	Type	Reusable O <sub>2</sub> Manifold
	Materials	Polysulfone Silicone Rubber
	Model #	SG 193-M
	Type	Reusable Adapter
	Size	22mm O.D. x 30mm I.D.
	Materials	Polysulfone

### EXPLANATION OF SYMBOLS

	Federal law restricts this device to sale by or on the order of a physician.		Non-Sterile
	Not manufactured with Natural Rubber Latex		Not manufactured with Di(2-ethylhexyl) phthalate (DEHP)
	Manufacturer		Single Patient Use

**Instrumentation**  
Industries, Inc.

2990 Industrial Blvd. • Bethel Park, PA 15102  
 US Toll Free: 1-800-633-8577  
 Business: 1-412-854-1133  
 US Toll Free Fax: 1-877-633-8661  
 Fax: 1-412-854-5668  
 E-mail: sales@iiimedical.com  
 www.iiimedical.com

Instrumentation  
Industries, Inc.

# O<sub>2</sub> Enrichment Kit

Installation  
& Usage  
Directions

**SG 066**



**Reusable**



## A Indications for Use

The SG 066 Oxygen (O<sub>2</sub>) Enrichment Kit is designed for use on certain ventilators to provide higher concentrations of oxygen to the patient by mixing room air with an oxygen source at the proximal airway.

## B Contraindications

None known.

## C Cautions & Warnings

**Warning:** Never operate the SG 066 O<sub>2</sub> Enrichment Kit without securing it to a ventilator.


**Warning:** The SG 066 O<sub>2</sub> Enrichment Kit is designed to operate with a hospital grade oxygen supply. Pressure should not exceed 10 psi at the O<sub>2</sub> Manifold Inlet. Oxygen flow to the SG 066 O<sub>2</sub> Enrichment Kit should not exceed 10 liters per minute.

**Caution:** Ensure that the Intake Port is correctly located, so that it cannot be blocked.

## D SG 066 O<sub>2</sub> Enrichment Kit Directions for Set-Up

### Directions for Set-up:

- Slide the SG 066 O<sub>2</sub> Enrichment Kit onto the mounting rail using the SG 066-3 O<sub>2</sub> Manifold Bracket. ① Ensure that the Ambient Air Intake Port is facing toward the front of the ventilator.
- Connect the SG 193-M Adapter to the KC 024 Tubing Outlet. ②
- Connect Small Bore Tubing (not included) to a medical grade, pressure compensated O<sub>2</sub> Flowmeter (not included). ③
- Connect the other end of the Small Bore Tubing to the O<sub>2</sub> Inlet at the top of the O<sub>2</sub> Manifold. ④
- Connect the KC 024 Tubing/SG 193-M Adapter Assembly to Inlet Filter on Ventilator. ⑤



## E SG 066 O<sub>2</sub> Enrichment Kit Directions for Operation

A physician must determine if there is a need for supplemental oxygen.

The SG 066 O<sub>2</sub> Enrichment Kit can deliver up to 100% oxygen concentration at the proximal airway.



- Determine the Minute Volume using the following formula:

$$\dot{V}_e = (\text{BPM})(V_t)$$

$\dot{V}_e$  = Minute Volume in liters/minute  
 $V_t$  = Tidal Volume  
 BPM = Breaths Per Minute

- Determine the Inspiratory:Expiratory (I:E) ratio (1:1 or 1:2).
- After determining the I:E ratio, use the tables below:
  - Locate the percentage of oxygen required on the Y-axis.
  - Locate the predetermined minute volume on the diagram.
  - Determine the flow rate of oxygen needed to achieve the desired oxygen concentration by locating the point on the X-axis directly below the established intersection.
- Analyze the delivered oxygen concentration for accuracy.

**Note:** These tables are not applicable if an oxygen concentrator is in use.



**Note:** During normal operation of the SG 066 O<sub>2</sub> Enrichment Kit, the reservoir bag will partially collapse with each cycle. This is due to the method at which the oxygen is blended with the room air.