



Instrumentation Industries, Inc.

Products for:

Anesthesia
Respiratory Care
Emergency Medicine

- Airway Management
- Emergency Ventilation
- O₂ Delivery
- Patient/Facility Accessories
- Pressure Monitoring
- Ventilator Circuits

Providing Quality Medical Products for over **58 years!**



Instrumentation Industries, Inc.

Since 1967

Our Story

We **breathe easy when the world breathes easier**...are watchwords for our small company. We saw the need for user-friendly products for medical/respiratory professionals and went to work to develop and produce practical solutions.

We manufactured just one product about **57 years ago** which we advertised on a postcard. Our latest catalog features over **200** products that contribute to the value, quality, comfort and safety of the healing process.

Our small and innovative production team continues to drive our mission to provide the very best respiratory/anesthesia products for patient care. Our employees are dedicated and most of them have worked with us for many years. Our motivation is knowing that our customers depend on us to produce reliable products.

Mission Statement

Instrumentation Industries, Inc. manufactures and distributes **respiratory** and **anesthesia devices** for the medical community. We provide practical solutions that complement respiratory and anesthesia applications.

Our focus is to contribute to the value, quality, comfort, and safety of the healing process. This shared sense of purpose encompasses all aspects of our organization.

Quality Policy

Instrumentation Industries, Inc. strives to provide our customers with safe, reliable and effective medical devices that meet or exceed their expectations. This is to be accomplished in accordance with all relevant national and international regulations.

We are committed to maintaining the effectiveness of the Quality Management System and the establishment, periodic review, and achievement of our Quality Objectives.

Instrumentation Industries, Inc. is certified to **ISO 13485:2016** and conforms to the requirements of **MDSAP** for the **USA** and **Canada**.



Environmental Policy

Instrumentation Industries, Inc. is committed to the protection of the environment and the health and safety of its employees, customers and neighbors. At III, our goals are continuous improvement, finding ways to do more with less, to reduce & recycle our waste, and to become efficient in every way possible.

All of our printed materials are printed in the USA on recycled paper with soy-based inks.

Visit iiimedical.com or call 1-800-633-8577 (inside USA) 1-412-854-1133 (outside USA) for our **REACH**, **RoHS**, and **Conflict Mineral Statements**.



Edward C. Horey, Jr.
President

5 ways to order:

1. Online: iiimedical.com

- ➔ Orders can be securely placed **online 24 hours a day**.
- ➔ Please include a purchase order number (*if available*) when ordering.

2. Phone: **1-800-633-8577** (*inside USA*) — **1-412-854-1133** (*outside USA*)

- ➔ Our Customer Service Representatives are available to take your orders **Monday through Friday, 8:30 A.M. to 4:00 P.M. EST.**
- ➔ Please have a purchase order number available when ordering.
- ➔ Volume discounts may be available. Please inquire at the time you place your order.

3. Fax: **1-877-633-8661** (*inside USA*) — **1-412-854-5668** (*outside USA*)

- ➔ Previous customers can fax orders **24 hours a day**.
- ➔ Please include a purchase order number and customer number when ordering.

4. Email: sales@iiimedical.com

- ➔ Email our friendly Customer Service 24 hours a day and we will get back to you to take your order.

5. EDI: ghx.com (PO850)

- ➔ We are a member of **GHX Healthcare Exchange**.
- ➔ Call **1-800-633-8577** (*inside USA*) — **1-412-854-1133** (*outside USA*) for information: (**850/PO810/Inv EDI Connection**)

You can also order from one of our Dealers.

- ➔ To find one of our worldwide dealers, visit iiimedical.com/dealers/.
or call **1-800-633-8577** (*inside USA*) — **1-412-854-1133** (*outside USA*)

Ordering Information:

- ➔ We accept **American Express, Visa or MasterCard.**



Shipping / Terms:

- ➔ **F.O.B. Bethel Park, PA, USA.** Prepaid and added to invoice.

Satisfaction Guarantee:

- ➔ All **Instrumentation Industries, Inc.** products are manufactured to exacting specifications, in accordance with all relevant national and international regulations, and are carefully inspected before shipment.
- ➔ Our products are fully guaranteed against defects in material and workmanship.
- ➔ Products will be replaced if found defective at the discretion of **Instrumentation Industries, Inc.**

Returns:

- ➔ Returns will not be accepted without an **Authorized Return Number**.
- ➔ Please have **Invoice Number** when calling to obtain the Authorized Return Number.
- ➔ Items returned will be subject to a **25% restocking fee**.
- ➔ Items must be returned within **30 days** in the original packaging.
- ➔ **Custom Orders** are not eligible for returns.

Technical Assistance:

- ➔ Our experienced staff can help with specific problems or questions regarding the use, function and application of our products. call **1-800-633-8577** (*inside USA*) — **1-412-854-1133** (*outside USA*)

General Product Information

All products are packaged as **Non-Sterile** (except ACU 1060.1 & ACU 1060.2).

Always follow the Sterilization and High Level Disinfection Recommendations included in the product's package insert, if applicable.

Visit iiimedical.com for:

- An up-to-date listing of **Sterilization and High Level Disinfection Recommendations**
- **Latex, DEHP, & BPA** information
- **Package Insert/User Instructions**, where applicable
- **Current pricing**

Product Symbols used in this catalog






	Reusable Products are intended for repeated use. Follow facility guidelines for reprocessing.
	Single Patient Use Products are intended to be used on a single patient only.
	Single Patient Reuse Products are intended to be used to repeatedly administer medications to one single patient while installed in a breathing circuit.
	Products with this symbol need a prescription to purchase. Call our customer service for information.
	Products with this symbol pose unacceptable risks to the patient, medical staff or other persons within the magnetic resonance imaging (MRI) environment.

Table of Contents

Product	Page
Adapter — Inline Aerosol	1
Adapters — Barbed Tubing (<i>"Christmas Tree" Connectors</i>)	25
Adapters — Swivel Barbed Tubing (<i>Swivel "Christmas Tree" Connectors</i>)	25
Adapter — Flowmeter/Air Compressor	27
Adapters — Gauge	18
Adapters — General	4
Adapters — MDI (<i>Metered Dose Inhaler</i>)	1
Adapters — Universal Inline MDI (<i>Metered Dose Inhaler</i>)	2
Adapters — Universal Inline MDI Kit (<i>Metered Dose Inhaler</i>)	2
Adapters — Inline MDI (<i>Metered Dose Inhaler</i>)	3
Adapters — Pressure	6
Adapters — Reducing	5
Adapters — Specialty	3
Adapters — Tee	4
Airways — Oropharyngeal	21
Anesthesia Flush Plug	30
Bag Hangers	27
Bags — Breathing	28
Brackets — Flexible Support Arms	12
Breast Marker	29
Breathing Bags	28
Bronchoscope Swivels & Caps	9
Cap — MDI Adapter	3
Cap — Patient Wye Connector Temp. Probe Port	8
Caps — Bronchoscope Swivel	9
Caps — Pressure Adapter	6
"Christmas Tree" Connectors — <i>Barbed Tubing Adapters</i>	25
Swivel "Christmas Tree" Connector — <i>Swivel Barbed Tubing Adapter</i>	25
Clamps — Gauge	18
Connector — 5-in-1	27
Connectors — 90° Elbow	6
Connectors — Hose	7
Connectors — Patient Wye	8
Connectors — Silicone Rubber	7
Cuff Pressure Monitors & Adapter	24
Endotracheal Tube Changers	23
Exhalation Valves	10
Exhalation Valves — Replacement Parts	10
Flexible Support Arms	12
Flexible Support Arms — Brackets	12
Flexible Support Arms — Accessories	13
Flowmeter/Air Compressor Adapter	27

Product	Page
Flowmeter Security Devices	29
Gauge — Pressure	16
Gauge Adapters	18
Gauge Clamps	18
Gauge Cover	19
Gauge Crystals	16, 17, 19
Gauge Stand	19
Gauges — Vacuum/Pressure	16
Hose Assemblies — High Pressure	28
Humidifier Guard	15
Knobs — Replacement	12
Laryngoscope Blades & Handles	22
Laryngoscope Replacement Bulbs	22
Manual Jet Ventilators & Replacement Parts	20
Mask/Filter Storage Boxes & Replacement Seals	30
Metered Dose Inhaler (MDI) Adapters	1-3
Metered Dose Inhaler Adapter Cap	3
NIF (Negative Inspiratory Force) Meters	17
Non-Rebreathing Tee Valve	9
Nose Clip	29
O₂ Instant Flow Valves	26-27
O₂ Instant Flow Valves Bag Hangers	27
O₂ Instant Flow Valve Wall Mount	27
O₂ Tank Wrench	29
One-Way Valves	11
Oropharyngeal Airways	21
Oxygen/No Smoking Signs	31
PEEP Valves	10
Pilot Tube Repair Kit	24
Safety Strap	13
Screwdriver — Reversible	17
SeCUREacord — Cord Covers	31
Silicone Rubber Connectors	7
Temp. Probe Port Cap — Patient Wye Connector	8
Test Lungs	11
Thermometer — LCD	9
Transtracheal Catheters	21
Vacuum/Pressure Gauges	16
Ventilator Tubing — Accessories	13
Ventilator Tubing — Cuffed & Cuffless	14-15
Water Traps	11

Inline Aerosol Adapter

Inline Aerosol Adapter allows a Combivent® Respimat® Inhaler to be adapted to a ventilator circuit. Immediately before the ventilator cycle, activate the Combivent® Respimat® Inhaler. The medication then becomes suspended in the Inline Aerosol Adapter until the inspiratory cycle of the ventilator carries it to the patient.

RTC 26-C

Inline Aerosol Adapter

- **Single Patient Reuse**
(Intended to be used to deliver Combivent® (ipratropium bromide and albuterol) Respimat® medication multiple times to a single patient while installed in a breathing circuit.)
- Designed to fit in the **Inspiratory Limb** of the patient circuit
- Combivent® Respimat® Inhaler may be placed on **RTC 26-C Inline Aerosol Adapter** in either orientation
- **Clear Body** allows viewing of suspended medication after activation
- **Arrows** clearly mark direction of airflow
- Tight-fitting **Thermoplastic Rubber Cap** maintains proper patient ventilation between treatments
- Polystyrene Butadiene
- Patent Pending
- Packaged: 25 or 50



Rx ONLY


SPR SINGLE PATIENT REUSE

MADE IN USA!

For use with
Combivent® Respimat® inhaler

Combivent® is a registered trademark of Boehringer Ingelheim Pharmaceuticals, Inc. Copyright © 2016 Boehringer Ingelheim International GmbH.

Respimat® is a registered trademark and used under license from Boehringer Ingelheim International GmbH.



[Click to View](#)
 — Installation & Usage Directions
 — Consignes d'installation et d'utilisation

[Buy Online](#)

MDI Adapters

MDI Adapters allow for delivery of aerosol medication from **Standard MDI Canisters**.

- Delivers medication directly via **Endotracheal Tube** or through **Fiber-Optic Bronchoscopy Adapters**
- Accepts canisters with **metal** or **plastic** tips

Rx ONLY

MADE IN USA!

For use with
Standard Canister

RTC 21-A MDI Adapter

- Polyetherimide
- Packaged: 5 or 10

R REUSABLE

[Buy Online](#)


RTC 21-D MDI Adapter

- Polystyrene Butadiene
- Packaged: 50 or 100

SPU SINGLE PATIENT USE

[Buy Online](#)

[Click to View](#)
 — Installation & Usage Directions



Inline MDI Adapters — pg3

Universal Inline MDI Adapters

Universal Inline MDI Adapters direct aerosol medication mainstream for ventilator patients.


- 2 types: for 22mm I.D. Ventilator Circuits & for 15mm I.D. Ventilator Circuits
- **Single Patient Reuse**
(Can be used to repeatedly administer medications to one single patient while installed in a breathing circuit.)
- Accepts MDI Canisters with a dose counter
- Advances certain dose counters
- Accepts Standard MDI Canisters with metal or plastic tips
- **Arrows** clearly mark direction of airflow
- Tight-fitting Thermoplastic Rubber Cap maintains proper patient ventilation between treatments
- Polystyrene Butadiene
- Packaged: 25 or 50



RTC 24-V
Adult Universal
Inline MDI Adapter

[Buy Online](#)

[Click to View](#)
— Installation & Usage Directions
— Consignes d'installation et d'utilisation



RTC 24-VP
15mm Universal
Inline MDI Adapter



[Buy Online](#)

[Click to View](#)
— Installation & Usage Directions

How they work

RTC 24-V/VP
Universal Inline
MDI Adapter

Both **Standard MDI Canisters** and **MDI Canisters with a dose counter** can be activated by **RTC 24-V/VP**



Rx ONLY

MADE IN USA!


SPR
SINGLE PATIENT REUSE

Adult Universal Inline MDI Adapter Kit

Adult Universal Inline MDI Adapter Kit creates a reservoir for our RTC 24-V Adult Universal Inline MDI Adapter.

RTC 24-V KIT
Adult Universal Inline MDI Adapter Kit

- Single Patient Reuse
(Can be used to repeatedly administer medications to one single patient while installed in a breathing circuit.)
- Can be used at the endotracheal tube or the inspiratory limb of 22mm I.D. circuits.
- Packaged: 10
- Includes:
 - A RTC 24-V — Adult Universal Inline MDI Adapter — Polystyrene Butadiene w/Thermoplastic Rubber Cap
 - B KC 2124 — Adapter — Polystyrene Butadiene
 - C KC 004-5 — Tubing — Polyolefin Plastomer





[Buy Online](#)

[Click to View](#)
— Installation & Usage Directions
— Consignes d'installation et d'utilisation

How it works

RTC 24-V KIT
Adult Universal Inline MDI
Adapter Kit



Rx ONLY

SPR
SINGLE PATIENT REUSE


MADE IN USA!



Inline MDI Adapters

Inline MDI Adapters direct aerosol medication mainstream from **Standard MDI Canisters** for ventilator patients.

- Tight-fitting **Cap** maintains proper patient ventilation between treatments
- **Arrows** clearly mark direction of airflow
- **6 inches of tubing** can be added to create a chamber when applicable
- Accepts canisters with **plastic** or **metal** tips



Rx ONLY

MADE IN USA!



[Click to View](#)

- Installation & Usage Directions
- Consignes d'installation et d'utilisation

RTC 15-D
Pediatric Inline MDI Adapter

- Thermoplastic Rubber Cap
- Polystyrene Butadiene
- Packaged: 50 or 100

Buy Online



RTC 15-P
Pediatric Inline MDI Adapter

- Thermoplastic Rubber Cap
- Polystyrene Butadiene
- Packaged: 50 or 100

Buy Online



RTC 22-A
Adult Inline MDI Adapter

- Silicone Rubber Cap
- Polyetherimide
- Packaged: 5 or 10

Buy Online



RTC 22-D
Adult Inline MDI Adapter

- Thermoplastic Rubber Cap
- Polystyrene Butadiene
- Packaged: 50 or 100

Buy Online



RTC 23-A
Adult Inline MDI Adapter

- Silicone Rubber Cap
- Polyetherimide
- Packaged: 5 or 10

Buy Online



RTC 23-D
Adult Inline MDI Adapter

- Thermoplastic Rubber Cap
- Polystyrene Butadiene
- Packaged: 50 or 100

Buy Online



RTC 17-1
Inline MDI Adapter Replacement Cap

- For RTC 22-A & RTC 23-A
- Silicone Rubber



R
REUSABLE

MADE IN USA!

- Packaged: 5 or 50

Buy Online

Specialty Adapters

Specialty Adapters increase treatment delivery options in special circumstances.

MADE IN USA!

BE 105-2222
Adult Swivel Adapter

- For CPAP/BiPAP systems
- Polypropylene
- Packaged: 10 or 100

SPU
SINGLE PATIENT USE

Buy Online



BE 106
Adult Adapter w/Vent

- Polycarbonate
- Packaged: Individually

R
REUSABLE

Buy Online



BE 120-AKW
45° Adapter

- Polystyrene Butadiene
- Packaged: 100

SPU
SINGLE PATIENT USE

Buy Online



BE 120-22D
Mouthpiece Adapter

- Polystyrene Butadiene
- Packaged: 100

SPU
SINGLE PATIENT USE

Buy Online



General Adapters

General Adapters are the same size on both ends for circuit customization when needed.

MADE
IN
USA!

Adult

KC 105
Adult
General Adapter

- Polycarbonate
- Packaged: 1 or 10



Buy Online

KC 2124
Adult
General Adapter

- Polystyrene Butadiene
- Packaged: 25 or 100



Buy Online

BE 105-2X22
Adult
General Adapter

- Polycarbonate
- Packaged: 5 or 25



Buy Online

BE 105-2X22D
Adult
General Adapter

- Polystyrene Butadiene
- Packaged: 20 or 100



Buy Online

Note: Transparent light blue product inventory in limited quantity. When depleted, product will be clear translucent, as shown.

Neonatal & Pediatric

KC 1010
Neonatal
General Adapter

- Polystyrene Butadiene
- Packaged: 25 or 100



Buy Online

SG 1010
Neonatal
General Adapter

- Polysulfone
- Packaged: 1 or 10



Buy Online

KC 1515
Pediatric
General Adapter

- Polystyrene Butadiene
- Packaged: 25 or 100



Buy Online

We also have other Adapters not specifically for ventilator circuits!

Gauge Adapters page 18

Flowmeter/
Air Compressor Adapter page 27

Tee Adapters

Tee Adapters help add versatility to patient circuits. Available in seven types to allow for various configurations during patient circuit assembly.

MADE
IN
USA!

BE 125
Tee Adapter

- Polycarbonate
- Packaged: 1 or 10



Buy Online

BE 125-1
Tee Adapter

- Polycarbonate
- Packaged: Individually



Buy Online

BE 125-1T
Tee Adapter

- Polycarbonate
- Packaged: 1 or 10



Buy Online

BE 125-1TD
Tee Adapter

- High Density Polyethylene
- Packaged: 10, 20 or 100



Buy Online

BE 16-2
Tee Adapter

- Polycarbonate
- Packaged: Individually



Buy Online

BE 117-W
Tee Adapter

- High Density Polyethylene
- Packaged: 10



Buy Online

Note: Blue product inventory in limited quantity. When depleted, product will be translucent white.

Reducing Adapters

Reducing Adapters are available in various configurations. **MADE IN USA!**

BE 190
Reducing Adapter
 ■ Polycarbonate
 ■ Packaged: 1 or 10




[Buy Online](#)

BE 191
Reducing Adapter
 ■ Polycarbonate
 ■ Packaged: 1 or 10





[Buy Online](#)

BE 192
Reducing Adapter
 ■ Polycarbonate
 ■ Packaged: 1 or 10





[Buy Online](#)

BE 193
Reducing Adapter
 ■ Polycarbonate
 ■ Packaged: 1 or 10




[Buy Online](#)

BE 139-105
Reducing Adapter
 ■ Polycarbonate
 ■ Packaged: 1 or 10





[Buy Online](#)

BE 195-1522A
Reducing Adapter
 ■ Polyetherimide
 ■ Packaged: 1 or 10





[Buy Online](#)

KC 2119
Reducing Adapter
 ■ Polystyrene Butadiene
 ■ Packaged: 25





[Buy Online](#)

KC 2126
Reducing Adapter
 ■ Polystyrene Butadiene
 ■ Packaged: 25 or 100





[Buy Online](#)

NB 105-6
Reducing Adapter
 ■ Polycarbonate
 ■ Packaged: Individually





[Buy Online](#)

Adult to Pediatric or Neonatal

SG 101
Adult to Pediatric Reducing Adapter
 ■ Polysulfone
 ■ Packaged: 1 or 10





[Buy Online](#)

SG 101-D
Adult to Pediatric Reducing Adapter
 ■ Polystyrene Butadiene
 ■ Packaged: 25 or 100





[Buy Online](#)

SG 102
Adult to Neonatal Reducing Adapter
 ■ Polysulfone
 ■ Packaged: 1 or 10




[Buy Online](#)



Our product quality is rigorously monitored throughout the entire manufacturing process.

Exceptional Product Quality

Pressure Adapters

Pressure Adapters provide easy access for pressure measuring.

- ▶ For use in Ventilator, BIPAP, or CPAP circuits
- ▶ Can be placed in any Adult, or Pediatric circuit

BE 105-40
Pressure Adapter

- Capped Barbed Port
Fits 3/16" I.D. tubing
- Silicone Rubber Port Cap
- Polycarbonate
- Packaged: Individually



[Buy Online](#)



15mm I.D.
22mm O.D.

MADE
IN
USA!

BE 105-101-40
Pressure Adapter

- Capped Barbed Port
Fits 3/16" I.D. tubing
- Silicone Rubber Port Cap
- Silicone Rubber Connector
- Polycarbonate
- Packaged: Individually

[Buy Online](#)



22mm I.D.

BE 139
Pressure Adapter

- Barbed Port
Fits 3/16" I.D. tubing
- Polycarbonate
- Packaged: 1 or 10



[Buy Online](#)



15mm I.D.
22mm O.D.

Page
16

Attach one of our
Vacuum/Pressure Gauges
for easy inline pressure measuring!

KC 5
Pressure Adapter

- Luer Lock Port
- Polycarbonate
- Packaged: Individually



[Buy Online](#)



15mm I.D.
22mm O.D.

KC 6PA
Pressure Adapter

- Capped Luer Lock Port
- Polystyrene Butadiene
- Packaged: 25 or 100



[Buy Online](#)



15mm I.D.

MS 101-40
Pressure Adapter

- Barbed Port
Fits 3/16" I.D. tubing
- Polycarbonate
- Packaged: Individually



[Buy Online](#)



15mm I.D.

NP 130-22
Pressure Adapter

- Capped Barbed Port
Fits 1/4" I.D. tubing
- Silicone Rubber Port Cap
- Polycarbonate
- Packaged: Individually



[Buy Online](#)



22mm I.D.

Replacement Caps

- Silicone Rubber



MADE
IN
USA!

BE 40A
Replacement Cap

- Fits 1/4" - 9/32" barbs
- Replacement for: NP 130-22
- Packaged: 5 or 50

[Buy Online](#)



BE 40H
Replacement Cap

- Fits 3/16" - 7/32" barbs
- Replacement for: BE 105-101-40 & BE 105-40
- Packaged: 5 or 50

[Buy Online](#)



90° Elbow Connectors

MADE
IN
USA!

90° Elbow Connectors are available with a luer lock port or a barbed port to allow for gas sampling or pressure monitoring.

Adapters vs Connectors

A Adapters generally connect different sizes together.

C Connectors generally connect items of the same size.

KC 18901-K
90° Elbow Connector

- Polystyrene Butadiene
- Packaged: 50



[Buy Online](#)



15mm O.D.

KC 18901-L
90° Elbow Connector

- Capped Luer Lock Port
- Polycarbonate
- Packaged: 10 or 20



[Buy Online](#)



15mm O.D.

KC 18901-O
90° Elbow Connector

- Barbed Port
Fits 3/16" I.D. tubing
- Polycarbonate
- Packaged: 10 or 20



[Buy Online](#)



15mm O.D.

KC 18901-P
90° Elbow Connector

- Polycarbonate
- Packaged: 10 or 20



[Buy Online](#)



15mm O.D.

Silicone Rubber Connectors

Silicone Rubber Connectors help to join circuit fittings when they don't quite fit. Dissimilar sizes, or accessories of the same size, can be used together.

Make use of "odd" size accessories.

- Available in **Adult, Pediatric, & Neonatal**
- Remains pliable after sterilization
- Autoclavable
- Silicone Rubber
- Packaged: 10 or 50

MADE
IN
USA!

BE 101-00-1426

Tapered
Silicone Rubber
Connector

- Adult
- Large end also fits flanged 22mm O.D.

R
REUSABLE

Buy Online



27mm I.D.
22mm I.D.

BE 101-00-3443

Straight
Silicone Rubber
Connector

- Adult

R
REUSABLE

Buy Online



22mm I.D.
22mm I.D.

PF 001

Reducing
Silicone Rubber
Connector

- Adult
to Pediatric

R
REUSABLE

Buy Online



22mm I.D.
15mm I.D.

PF 1510

Reducing
Silicone Rubber
Connector

- Pediatric
to Neonatal

R
REUSABLE

Buy Online



15mm I.D.
10mm I.D.

PF 1515

Straight
Silicone Rubber
Connector

- Pediatric

R
REUSABLE

Buy Online



15mm I.D.
15mm I.D.

PF 1010

Straight
Silicone Rubber
Connector

- Neonatal

R
REUSABLE

Buy Online



10mm I.D.
10mm I.D.

Hose Connectors

1/4" Tapered Barbed Connection

Hose Connectors change a 15mm, 22mm or 30mm I.D. low pressure gas outlet source to a 1/4" I.D. small bore tubing connection. Also available, a 15mm O.D. x 1/2" hose barb connector (ADP 150) that fits 10mm I.D. tubing.

MADE
IN
USA!

BE 145-1522P

Hose Connector

- Polycarbonate
- Packaged:
10 or 50

R
REUSABLE

Buy Online



15mm I.D.
22mm O.D.

KC 145-1522D

Hose Connector

- Polystyrene
Butadiene
- Packaged:
25 or 100

SPU
SINGLE
PATIENT
USE

Buy Online



15mm I.D.
22mm O.D.

BE 145-2P

Hose Connector

- Polycarbonate
- Packaged:
10 or 50

R
REUSABLE

Buy Online



15mm O.D.

BE 145-2230P

Hose Connector

- Polycarbonate
- Packaged:
10 or 50

R
REUSABLE

Buy Online



22mm I.D.
30mm O.D.

KC 145-2230D

Hose Connector

- Polystyrene
Butadiene
- Packaged:
25 or 100

SPU
SINGLE
PATIENT
USE

Buy Online



22mm I.D.
30mm O.D.

ADP 150

Hose Connector

- Fits cuffed
10mm I.D.
Tubing
- Polycarbonate
- Packaged:
10 or 50

R
REUSABLE

Buy Online



1/2" O.D.

1/2" O.D.
Hose Barb

Ventilator Tubing

In various
lengths of
adult,
pediatric,
& neonatal
sizes



pages 14-15

Patient Wye Connectors

Reusable Patient Wye Connectors are available with pressure ports and/or temperature probe ports for easy monitoring of the circuit at the patient wye. **MADE IN USA!**

BE 127
Patient Wye Connector

- No monitoring ports
- Polycarbonate
- Packaged: 1 or 10



[Buy Online](#)



BE 127-40
Patient Wye Connector

- Pressure Port Fits 1/4" I.D. tubing
- Polycarbonate
- Packaged: 1 or 10



[Buy Online](#)




BE 127-7.5
Patient Wye Connector

- 7.5mm I.D. Temperature Probe Port
- Polycarbonate
- Packaged: 1 or 10



[Buy Online](#)




BE 127-7.5-40
Patient Wye Connector

- 7.5mm I.D. Temperature Probe Port
- Pressure Port Fits 1/4" I.D. tubing
- Polycarbonate
- Packaged: 1 or 10



[Buy Online](#)



Temperature Probe Port Cap

Temperature Probe Port Cap attached to a wye connector helps maintain proper ventilation.

- Fits: BE 127-7.5 & BE 127-7.5-40 Temperature Probe Ports
- Thermoplastic Rubber
- Packaged: 5 or 50



[Buy Online](#)



CAP-S

Single Patient Use Patient Wye Connectors are available with monitoring ports at the patient wye. **MADE IN USA!**


Note: Transparent light blue product inventory in limited quantity. When depleted, product will be clear translucent.

KC 211-D
60° Patient Wye Connector

- 60° Wye
- 18mm I.D. 22mm O.D. Port
- Polystyrene Butadiene
- Packaged: 50 or 100



[Buy Online](#)




KC 212-D
90° Patient Wye Connector

- 90° Wye
- 18mm I.D. 22mm O.D. Port
- Polystyrene Butadiene
- Packaged: 50 or 100



[Buy Online](#)



KC 213-D
60° Patient Wye Connector

- 60° Wye
- No Monitoring Ports
- Polystyrene Butadiene
- Packaged: 50 or 100



[Buy Online](#)



Patient circuits need

3

things...

1. Tubing
2. Support
3. Accessories

We provide all three!

Bronchoscope Swivels

Bronchoscope Swivels allow access to the patient's airway, while maintaining ventilation.

- Two **Interchangeable Silicone Rubber Caps**
 - One cap allows access for a **Fiber-Optic Bronchoscope**
 - One cap allows access for a **Suction Catheter**
- Replacement Caps available
- Polyetherimide
- Packaged: 1 or 10



BE 105-7
Bronchoscope Swivel
15mm O.D. Equipment Connection
15mm I.D. Patient Connection

[Buy Online](#)



BE 105-8
Bronchoscope Swivel
15mm O.D. Equipment Connection
15mm I.D. 22mm O.D. Patient Connection

[Buy Online](#)

Replacement Caps

- Silicone Rubber



BE 105-7-1
Fiber-Optic Cap

- Packaged: 5



[Buy Online](#)



RTC 17-1
Suctioning Cap

- Packaged: 5 or 50



[Buy Online](#)

Bronchoscope Double Swivel

Bronchoscope Double Swivel has a **single cap** for both **fiber-optic** and **suction** access to a patient's airway.

DS 7-45
Bronchoscope Double Swivel

- Both Ends Swivel
- Polypropylene
- Silicone Rubber Cap
- Packaged: 25 or 50



15mm O.D. Equipment Connection
15mm I.D. 22mm O.D. Patient Connection

[Buy Online](#)

LCD Thermometer

LCD Thermometer helps monitor the temperature of ventilator or heated aerosol circuits.

Note: LCD readings are approximate and not a substitute for definitive temperature monitoring.

KM 030
LCD Thermometer

- LCD displays 28°, 31°, 34°, 37°, & 40° Celsius markings.
- High Density Polyethylene
- Packaged: 50



22mm I.D.


[Buy Online](#)

Non-Rebreathing Tee Valve

Non-Rebreathing Tee Valve isolates airflow and **decreases** the possibility of **nosocomial infection**. It can also reduce **respirometer** replacement costs. Low cracking pressure minimizes resistance to flow!

BE 117
Non-Rebreathing Tee Valve

- Polycarbonate
- Silicone Rubber Leaf Valves
- Packaged: 1 or 10




22mm O.D. Equipment Connection
15mm I.D. 22mm O.D. Patient Connection
15mm O.D. 22mm O.D. Oxygen/Air Connection

[Click to View](#)
Installation & Usage Directions

[Buy Online](#)

How it works

BE 117 — Non-Rebreathing Tee Valve



Magnetic PEEP Valves

Magnetic PEEP Valves are used with ventilators, resuscitators, or any system where PEEP or CPAP is required.

Caution:

- These are uni-directional PEEP Valves.
- Flow direction must be consistent with the expiratory flow.
- PEEP pressure may vary according to position.
- PEEP pressure must be monitored during use.

- Adjustable from 3 to 20 cm H₂O pressure
- Clearly marked flow direction arrows
- Transparent
- Polycarbonate, Polypropylene, Silicone Rubber, & Stainless Steel
- Packaged: Individually



BE 142
Magnetic PEEP Valve
w/Silicone Rubber Connector

[Buy Online](#)


BE 30-142
Magnetic PEEP Valve

[Buy Online](#)


How they work

BE 142/BE 30-142 — Magnetic PEEP Valve

When used with a Ventilator




When used with a Bag/Mask



Exhalation Valves

Exhalation Valves are available with various patient/inlet configurations. Replacement parts are available.

- Available in Regular or Low Resistance
- Arrows clearly mark direction of airflow
- Polyetherimide ■ Packaged: Individually



BE 30-115-B
Exhalation Valve
BE 30-115-BL
Low Resistance
Exhalation Valve

[Buy Online](#)

NS 30-115-U
Exhalation Valve
NS 30-115-UL
Low Resistance
Exhalation Valve

[Buy Online](#)


NS 30-115-5
Exhalation Valve
NS 30-115-5L
Low Resistance
Exhalation Valve

[Buy Online](#)


How they work

Exhalation Valves


BE 30-115-B
BE 30-115-BL
Side View



NS 30-115-5
NS 30-115-5L
Top View



NS 30-115-U
NS 30-115-UL
Side View



Exhalation Valve Replacement Parts

- Fits all of our Exhalation Valves

BE 215-2A — Diaphragm Assembly

- Silicone Rubber Diaphragm
- Polypropylene Stem & Ring
- Packaged: 1 or 10

[Buy Online](#)

BE 115-2N — Diaphragm Nut

- Nylon
- Packaged: 5 or 25

[Buy Online](#)



[Click to View](#)
Installation
& Usage
Directions

One-Way Valves


One-Way Valves are available in different configurations to allow various placement options in a patient circuit. Clearly marked **flow-direction arrows** help ensure proper placement in the circuit.

- Low Resistance
- **Arrows** clearly mark direction of airflow
- Silicone Rubber Leaf Valve
- Packaged: 1 or 10



Click to View

- Installation & Usage Directions
- Consignes d'installation et d'utilisation



BE 130
One-Way Valve

- Safety Lugs
- Polycarbonate

[Buy Online](#)



BE 130-22B
One-Way Valve

- Polysulfone

[Buy Online](#)



BE 131
One-Way Valve

- Safety Lugs
- Silicone Rubber Connector
- Polycarbonate

[Buy Online](#)



BE 130-23BB
One-Way Valve

- Polysulfone

[Buy Online](#)



Water Traps

Water Traps allow excess moisture from humidity to be collected. A spring-loaded plunger maintains circuit integrity while bottle is emptied.

- Available in 2 capacities: **4 oz. (125 ml)** & **8 oz. (250 ml)**
- Internal Spring: Steel
- Connector: Polyetherimide
- Bottle & Bottle Insert: Polypropylene
- Packaged: **Individually**



WT 125
Water Trap

- 4 oz. (125 ml)

[Buy Online](#)



WT 250
Water Trap

- 8 oz. (250 ml)

[Buy Online](#)



Replacement Parts

- Polypropylene
- Packaged: **Individually**



WT 125-BOT
4 oz. (125 ml)
Replacement Bottle

- For: WT 125

[Buy Online](#)



WT 125-BI
Replacement
Bottle Insert

- For: WT 125


[Buy Online](#)



WT 250-BOT
8 oz. (250 ml)
Replacement Bottle

- For: WT 250

[Buy Online](#)



WT 250-BI
Replacement
Bottle Insert

- For: WT 250

[Buy Online](#)

Test Lungs

Test Lungs are perfect for testing ventilators, resuscitators, etc. Not made with natural rubber latex.

- Available in 2 capacities: **750 ml** and **1-liter**
- Not for patient use
- Cannot be cleaned with alcohol
- Green Neoprene
- Packaged: 1 or 5



BE 150-5G
Test Lung

- 1-liter capacity

[Buy Online](#)



BE 150-5LF
Test Lung

- 750-ml capacity

[Buy Online](#)



Flexible Support Arms

Flexible Support Arms fit most ventilators or can be adapted to others by using one of our **Pole or Rail Mounts**.

- Available in 2 sizes, 24" & 30"
- Rotates 360°
- Sleeve Base locks in place
- Vinyl-coated segmented arm
- Enhances patient comfort & safety
- Can be adapted to most ventilators by using our **Flexible Support Arm Brackets**
- Anodized Aluminum & Steel
- Packaged: **Individually**


R REUSABLE **MADE IN USA!** [Click to View Installation & Usage Directions](#)

BE 122-24
Flexible Support Arm
▪ 24" & Extension for 38" total length

[Buy Online](#)

BE 122-30
Flexible Support Arm
▪ 30" & Extension for 44" total length

[Buy Online](#)



Replacement Parts

BE 122-S
Replacement Sleeve Base
▪ Fits: BE 122-24 & BE 122-30
▪ Anodized Aluminum
▪ Packaged: **Individually**

R REUSABLE **MADE IN USA!** [Buy Online](#)

KC 1
Replacement Spring Ball Clamp
▪ Fits: BE 122-24 & BE 122-30
▪ Vinyl-Coated Steel
▪ Packaged: **Individually**

R REUSABLE [Buy Online](#)

Flexible Support Arm Rail-Mount Brackets

Flexible Support Arm Rail-Mount Brackets add versatility by allowing the mounting of **Flexible Support Arms** to **Rails**.

- Anodized Aluminum & Steel
- Packaged: **Individually**

R REUSABLE **MADE IN USA!**



BE 114-5
Universal Rail-Mounted Bracket
▪ Permits rail mounting of **Flexible Support Arms**
▪ Does not fit on rails mounted flush to wall

[Buy Online](#)



RB 002
Compact Rail-Mounted Bracket
▪ Permits rail mounting of **Flexible Support Arms**
▪ For use with Fairfield-style rails only

[Buy Online](#)

Replacement Knobs

For Flexible Support Arms & O₂ Instant Flow Valves (pgs. 26 & 27)

- Steel
- Packaged: **Individually**

R REUSABLE **MADE IN USA!**

RB KNOB
Rail-Mounted Bracket Replacement Knob
▪ Fits: **RB 002**

[Buy Online](#)

BE 114-K
Replacement Knob
▪ Fits: **BE 3000-PA** (pg. 26)
BE 3000-PS (pg. 26)

[Buy Online](#)

BE 114-5K
Replacement Knob
▪ Fits: **BE 3000-RA** (pg. 27)
BE 3000-RS (pg. 27)
BE 114-5

[Buy Online](#)



Flexible Support Arm Tubing Clamps

Flexible Support Arm Tubing Clamps hold large bore tubing in place. Specially designed for our **Flexible Support Arms**.

- Two types are available:
 - **Single Clamp** in 1 size
 - **Double Clamp** in 2 sizes, 4" or 8"
- Phenolic Ball Knob
- Polypropylene
- Anodized Aluminum (*Double Tubing Clamps only*)
- Steel
- Packaged: **Individually**



Single Tubing Clamp



BE 114-5TB
Single Tubing Clamp

[Buy Online](#)

Double Tubing Clamps



BE 122-2TSB
4" Double Tubing Clamp

[Buy Online](#)



BE 122-2TB
8" Double Tubing Clamp

[Buy Online](#)

Ventilator Tubing Clamps

Swivel Tubing Clamp

Swivel Tubing Clamp attaches large bore tubing to bedrails.

BE 102
Swivel Tubing Clamp

- Fits tubing and rails 1/2" to 1" in diameter
- Swivels 360°
- Polypropylene
- Packaged: 20, 50, or 100



MADE IN USA!

[Buy Online](#)

Tubing Separator

Tubing Separator keeps small bore tubing fastened to, but separate from, large bore tubing.

KC 10
Tubing Separator

- Silicone Rubber
- Packaged: 10 or 100



MADE IN USA!

[Buy Online](#)

Safety Strap

Safety Strap allows any adapter, connector, or valve with **safety lugs** to be safely attached to a ventilator circuit. Secures ventilator circuit attachment to a patient's airway. Helps prevent accidental disconnects. **Not made with natural rubber latex.**

NP 105-1
Safety Strap

- 5 7/8" long
- 3/16" diameter center hole
- Neoprene
- Packaged: 10 or 50



MADE IN USA!

[Buy Online](#)

Do you have an idea for a new medical product?
We're always on the lookout for ideas that can help in patient treatment.

Contact us 1-800-633-8577 sales@iiimedical.com



Cuffed Ventilator Tubing



Cuffed Ventilator Tubing is transparent, non-conductive, and has a **corrugated spiral design** to prevent kinking. It has a **smooth internal diameter** which promotes laminar gas flow, and resists the trapping of water. Our ventilator tubing has industry standard connection diameters of **10, 15, or 22mm I.D.**




- No minimum on Stock Lengths!
- Packaged: **Individually**

[Click to View](#)

- 👉 — Installation & Usage Directions
- 👉 — Consignes d'installation et d'utilisation

Cuffed Ventilator Tubing

 <p>KC Series Adult Ventilator Tubing</p> <ul style="list-style-type: none"> ■ 22mm I.D. ■ Polyester Elastomer ■ Silicone Rubber Cuffs <p>Stock Lengths</p> <p>KC 006 - 6" KC 036 -36" KC 012 -12" KC 048 -48" KC 018 -18" KC 060 -60" KC 024 -24" KC 072 -72" KC 030 -30"</p> <p style="text-align: right;">Buy Online 👉</p>	 <p>KC-P Series Pediatric Ventilator Tubing</p> <ul style="list-style-type: none"> ■ 15mm I.D. ■ Polyester Elastomer ■ Silicone Rubber Cuffs <p>Stock Lengths</p> <p>KC 012-P -12" KC 048-P -48" KC 018-P -18" KC 060-P -60" KC 024-P -24" KC 072-P -72"</p> <p style="text-align: right;">Buy Online 👉</p>	 <p>KC-N Series Neonatal Ventilator Tubing</p> <ul style="list-style-type: none"> ■ 10mm I.D. ■ Polyester Elastomer ■ Silicone Rubber Cuffs <p>Stock Lengths</p> <p>KC 006-N - 6" KC 048-N -48" KC 018-N -18" KC 060-N -60" KC 024-N -24"</p> <p style="text-align: right;">Buy Online 👉</p>
---	---	--

Custom lengths of any style of **Ventilator Tubing** are available.
*Minimum order requirements apply.
 Call 1-800-633-8577 for info.*

All Patient Circuits need 3 things...

1. Tubing
2. Support
3. Accessories

Instrumentation Industries, Inc. offers all 3!
Always the highest quality products!

Cuffless Ventilator Tubing




<p>KC-3 Series Adult Ventilator Tubing</p> <ul style="list-style-type: none"> ■ 22mm I.D. ■ Polyester Elastomer <p style="text-align: right;">Buy Online 👉</p>	<p>Stock Lengths</p> <p>KC 006-3 - 6" KC 060-3 -60" KC 018-3 -18" KC 072-3 -72" KC 036-3 -36"</p>	
---	--	--

Clear Cuffless Ventilator Tubing

Clear Cuffless Ventilator Tubing is transparent, non-conductive and has a corrugated spiral design to prevent kinking. It has a smooth internal diameter which promotes laminar gas flow, and resists the trapping of water. Our ventilator tubing has industry standard connection diameters of 10, 15, or 22mm I.D.

- No minimum on Stock Lengths!
- Packaged: **Individually**

[Click to View](#)

- Installation & Usage Directions
- Consignes d'installation et d'utilisation


Clear Cuffless Ventilator Tubing		
		
<p>KC-5 Series Adult Ventilator Tubing</p> <ul style="list-style-type: none"> ■ 22mm I.D. ■ Polyolefin Plastomer 	<p>KC-5P Series Pediatric Ventilator Tubing</p> <ul style="list-style-type: none"> ■ 15mm I.D. ■ Polyolefin Plastomer 	<p>KC-5N Series Neonatal Ventilator Tubing</p> <ul style="list-style-type: none"> ■ 10mm I.D. ■ Polyolefin Plastomer
<p>Stock Lengths</p> <p>KC 004-5 - 4" KC 024-5 - 24" KC 006-5 - 6" KC 060-5 - 60" KC 018-5 - 18" KC 072-5 - 72"</p>	<p>Stock Lengths</p> <p>KC 006-5P - 6" KC 048-5P - 48" KC 012-5P - 12" KC 060-5P - 60" KC 024-5P - 24"</p>	<p>Stock Lengths</p> <p>KC 004-5N - 4" KC 006-5N - 6" KC 024-5N - 24"</p>
Buy Online	Buy Online	Buy Online



Humidifier Guard

NEW!

Our **New HG 100 Humidifier Guard** helps to protect humidifiers by allowing unwanted fluids to be collected and then easily and safely removed by suction. It installs onto the humidifier's outlet port. Accumulated condensation and patient secretions in the **inspiratory limb** of a ventilator circuit are easily and safely removed by activating the attached closed **Suction System**. A **Temperature Probe** can be inserted directly into the **HG 100 Humidifier Guard** to monitor circuit temperature.



- No potential pressure loss in ventilator circuit
- Specifically designed for use with humidifiers in ventilator circuits
- Clearly marked **Ports, Direction Arrows, and FULL Line**
- Maximum amount of collected fluids is **40ml**
- The HG 100 Humidifier Guard **CANNOT** be used during **Magnetic Resonance Imaging (MRI)** as some parts of the assembly are magnetic
- Packaged: 5 or 10 ■ Patent Pending

HG 100
Humidifier Guard

[Click to View](#)

- Installation & Usage Directions



[Buy Online](#)

Vacuum/Pressure Gauges

Vacuum/Pressure Gauges provide accurate vacuum or pressure monitoring.

- Accurate
- Sturdy
- Reliable
- Incremental graduations as low as 1cm H₂O
- Easily re-zeroed for accuracy
- Internal Components: Aluminum, Brass, Copper Alloy & Steel
- Case: Acrylonitrile Butadiene Styrene (ABS)
- Crystal: Acrylic
- Packaged: Individually



[Click to View](#)
 Installation & Usage Directions
 Consignes d'installation et d'utilisation



***NS-S Series Gauge info for California customers:**
 WARNING: These products can expose you to chemicals including Lead, which is known to the State of California to cause cancer and birth defects or other reproductive harm. For more information go to www.P65Warnings.ca.gov.

Replacement Crystal for Vacuum/Pressure Gauges

Use as replacement crystal for any Instrumentation Industries, Inc. NS Series Vacuum/Pressure Gauge.

NS-SC

- Acrylic
- Packaged: Individually



[Buy Online](#)



NS 30-PBS Vacuum/Pressure Gauge

- Reads ± 30cm H₂O in 1cm graduations
- Bottom 1/8" small bore tubing connection

[Buy Online](#)



NS 30-TRS Vacuum/Pressure Gauge

- Reads ± 30cm H₂O in 1cm graduations
- Rear 1/4" NPT connection

[Buy Online](#)



NS 60-PBS Vacuum/Pressure Gauge

- Reads ± 60cm H₂O in 1cm graduations
- Bottom 1/8" small bore tubing connection

[Buy Online](#)



NS 60-TBS Vacuum/Pressure Gauge

- Reads ± 60cm H₂O in 2cm graduations
- Bottom 1/4" NPT connection

[Buy Online](#)



NS 60-TRS Vacuum/Pressure Gauge

- Reads ± 60cm H₂O in 2cm graduations
- Rear 1/4" NPT connection

[Buy Online](#)




NS 120-TRS Vacuum/Pressure Gauge

- Reads ± 120cm H₂O in 4cm graduations
- Rear 1/4" NPT connection

[Buy Online](#)

How they work

NS -S Series Vacuum/Pressure Gauges



Pressure Gauge

Pressure Gauge provides accurate 0 to 120cm H₂O pressure monitoring.

NS 120P-TRS Pressure Gauge

- Reads 0 to 120cm H₂O in 2cm graduations
- Rear 1/4" NPT connection
- Internal Components: Aluminum, Brass, Copper Alloy & Steel
- Case: Acrylonitrile Butadiene Styrene (ABS)
- Crystal: Acrylic
- Packaged: Individually



[Buy Online](#)

**California customers see above info for NS 120P-TRS Pressure Gauge*

[Click to View](#)

Installation & Usage Directions
 Consignes d'installation et d'utilisation

NIF Meters (Negative Inspiratory Force)

NIF Meters (Negative Inspiratory Force) capture a patient's best effort with a resettable pointer.

- Accurate
- Sturdy
- Reliable
- Resettable Recording Needle
- Incremental graduations as low as 1cm H₂O
- Easily re-zeroed for accuracy
- Internal Components: Aluminum, Brass, Copper Alloy & Steel
- Case: Acrylonitrile Butadiene Styrene (ABS)
- Crystal: Acrylic
- Packaged: **Individually**



[Click to View](#)
Installation
& Usage Directions



NS-R Series NIF Meter info for California customers:
WARNING: These products can expose you to chemicals including Lead, which is known to the State of California to cause cancer and birth defects or other reproductive harm. For more information go to www.P65Warnings.ca.gov.



- NS 30-PBR**
NIF Meter
- Reads $\pm 30\text{cm H}_2\text{O}$ in 1cm graduations
 - Bottom 1/8" small bore tubing connection
 - Includes: BE 148-6 Gauge Adapter

[Buy Online](#)



- NS 30-TRR**
NIF Meter
- Reads $\pm 30\text{cm H}_2\text{O}$ in 1cm graduations
 - Rear 1/4" NPT connection
 - Includes: BE 149-1 Gauge Adapter

[Buy Online](#)



- NS 60-PBR**
NIF Meter
- Reads $\pm 60\text{cm H}_2\text{O}$ in 2cm graduations
 - Bottom 1/8" small bore tubing connection
 - Includes: BE 148-6 Gauge Adapter

[Buy Online](#)



- NS 60-TBR**
NIF Meter
- Reads $\pm 60\text{cm H}_2\text{O}$ in 2cm graduations
 - Bottom 1/4" NPT connection
 - Includes: BE 149-1 Gauge Adapter

[Buy Online](#)



- NS 60-TRR**
NIF Meter
- Reads $\pm 60\text{cm H}_2\text{O}$ in 2cm graduations
 - Rear 1/4" NPT connection
 - Includes: BE 149-1 Gauge Adapter

[Buy Online](#)



- NS 120-TRR**
NIF Meter
- Reads $\pm 120\text{cm H}_2\text{O}$ in 4cm graduations
 - Rear 1/4" NPT connection
 - Includes: BE 149-1 Gauge Adapter

[Buy Online](#)

Replacement Crystal for NIF Meters

Use as replacement crystal for any Instrumentation Industries, Inc. NS Series NIF Meter.

NS-RC

- Resettable Recording Needle
- Acrylic
- Packaged: **Individually**



[Buy Online](#)

Reversible Screwdriver

Reversible Screwdriver has one flathead end and one Phillips end...very handy for re-zeroing gauges or making fine adjustments.

EH 2

Reversible Screwdriver

- 5" long
- Pocket clip
- Plated Steel
- Packaged: 1 or 5




[Buy Online](#)

How they work

NS -R Series NIF Meters


1/8" Small Bore Tubing Connection

Small Bore Tubing (not included)
BE 148-6 Gauge Adapter (included)
Endotracheal Tube




1/4" NPT Connection

BE 149-1 Gauge Adapter (included)
Endotracheal Tube



1/4" NPT Connection

BE 148-1 Gauge Adapter (not included) pg. 22
Small Bore Tubing (not included)
BE 148-6 Gauge Adapter (not included) pg. 22
Endotracheal Tube



Gauge Adapters

Gauge Adapters allow for easy remote inspiratory measurement. Designed for all Instrumentation Industries, Inc. Vacuum/Pressure Gauges and NIF Meters.

- Polycarbonate
- Gauge not included
- Packaged: Individually



BE 148-1
Gauge Adapter
 ■ Changes a 1/4" NPT connection to a 1/8" or 3/16" small bore tubing connection

[Buy Online](#)

BE 148-6
Gauge Adapter
 ■ 15mm I.D./22mm O.D. x 1/8" or 3/16" small bore tubing connection for remote inspiratory monitoring

[Buy Online](#)

BE 149-1
Gauge Adapter
 ■ 15mm I.D./22mm O.D. x 1/4" NPT connection for inspiratory monitoring

[Buy Online](#)

Gauge Clamps

Gauge Clamps offer the convenience of mounting Vacuum/Pressure Gauges and NIF Meters in a variety of places.

- Mounts Rear 1/4" NPT-threaded gauges
- Small Bore Tubing Connection
- Precision-machined
- Anodized Aluminum & Steel
- Plated Brass Fittings
- Gauge not included
- Packaged: Individually



[Click to View](#)
Installation & Usage Directions

WARNING: These products can expose you to chemicals including Lead, which is known to the State of California to cause cancer and birth defects or other reproductive harm. For more information go to www.P65Warnings.ca.gov.

Gauge Pole Clamps



BE 148-5-1
Gauge Pole Clamp
 ■ Mounts gauges to poles with a diameter of 3/4"
 ■ See pg. 12 for Replacement Knob

[Buy Online](#)

BE 148-5-2
Gauge Pole Clamp
 ■ Mounts gauges to poles with a diameter of 3/4" to 1 1/2"
 ■ See pg. 12 for Replacement Knob


[Buy Online](#)

Gauge Clamp



BE 148-5-2.5
Gauge Clamp
 ■ Clamp opens 1" to 2 1/2"
 ■ Mounts gauges securely to poles, incubators, ventilators, etc.

[Buy Online](#)



Did you know?
 Our shipping materials are all made from recycled materials. All of our printed pieces are printed on recycled paper with soy inks.


Gauge Cover

Gauge Cover is designed to fit all of our **NS Series NIF Meters** and **Vacuum/Pressure Gauges** with the highest level of protection from falls, drops, vibrations, and other shocks.

NS-GC

Gauge Cover

- Unique rubber "Clam Shell" design
- Ribbed for maximum shock protection
- Melt-Processible Rubber
- Gauge not included
- Packaged: **Individually**



Gauge Stand

Gauge Stand keeps a gauge upright, stable and easy to view. Specially designed for use with our **NS Series -PBR, -PBS, -TRR, & -TRS Vacuum/ Pressure Gauges** and **NIF Meters**.

NS-GS

Gauge Stand

- "No-skid" rubber feet
- Gauges are easily inserted
- Polystyrene Butadiene
- Gauge not included
- Packaged: **Individually**



Accessories for Discontinued, Instrumentation Industries, Inc. Gauges with Metal Cases

Call 1-800-633-8577
or 1-412-854-1133 *Outside USA*

BE 147-GC
Blue Gauge Cover
w/Long Neck



- Fits Vacuum/Pressure Gauges:
BE 147, BE 147-30, BE 148, BE 148-30,
& BE 148-120
- Fits NIF Meters:
BE 147-149, BE 149, & BE 149-120
- Packaged: **Individually**

BE 147-GCV
Blue Gauge Cover
w/Short Neck



- Fits Vacuum/Pressure Gauges:
BE 147, BE 147-30, BE 148, BE 148-30,
& BE 148-120
- Fits NIF Meters:
BE 147-149, BE 149, & BE 149-120
- Packaged: **Individually**

BE 148-2
Screw-on Smooth Crystal



- Replacement Crystal for
Vacuum/Pressure Gauges:
BE 147, BE 147-30, BE 148, BE 148-30,
& BE 148-120
- Packaged: **Individually**

How are we doing?

Our #1 objective as a company is to provide medical professionals high-quality products that will help provide the best care.
feedback@iiimedical.com



Manual Jet Ventilators

Manual Jet Ventilators are lightweight, portable, durable, and essential for EMS, Anesthesia, and Emergency departments.

- Specially designed for use with our Transtracheal Catheters pg. 21
- Plated Brass Fittings
- Anodized Aluminum & Steel
- Packaged: **Individually**

BE 183-SUR

Manual Jet Ventilator w/Regulator & Gauge

- Reusable
 - Adjustable Regulator & Gauge*
 - On/Off Valve
 - 6-ft. High Pressure Hose Assembly
- Single Patient Use
 - 4-ft. Small Bore Tubing Assembly w/inline filter (see BE 183-2 below for replacement)

Caution: BE 183-SUR is designed to operate with pressure that does not exceed 50psi. Use Regulator & Gauge to assure operating pressure does not exceed 50psi.



Buy Online

[Click to View](#) Installation & Usage Directions
Consignes d'installation et d'utilisation



*BE 183-G Gauge info for California customers:
WARNING: This product can expose you to chemicals including Lead, which is known to the State of California to cause cancer and birth defects or other reproductive harm. For more information go to www.P65Warnings.ca.gov.

BE 183-SU

Manual Jet Ventilator

- Reusable
 - On/Off Valve
 - 6-ft. High Pressure Hose Assembly
- Single Patient Use
 - 4-ft. Small Bore Tubing Assembly w/inline filter (see BE 183-2 below for replacement)

Caution: BE 183-SU is designed to operate directly from a gas source that does not exceed 50psi.



Buy Online

[Click to View](#) Installation & Usage Directions
Consignes d'installation et d'utilisation

Replacement Parts for Manual Jet Ventilators



BE 183-2
4-ft. Small Bore Tubing Assembly w/ Inline Filter

- Medical Grade Tygon® Tubing
- For BE 183-SU & BE 183-SUR
- Packaged: **Individually**



Buy Online



BE 183-6
6-ft. High Pressure Hose Assembly

- For BE 183-SU & BE 183-SUR
- Standard Plated Brass Female DISS fittings on both ends
- Medical Grade PVC Tubing
- Packaged: **Individually**



MADE IN USA!

Buy Online



BE 183-R
Regulator & Gauge

- For BE 183-SUR
- Packaged: **Individually**



Buy Online



BE 183-G
Gauge

- For BE 183-SUR
- 0-60 psi
- Must be used with regulator
- Packaged: **Individually**



Buy Online

*California customers see above info for BE 183-G Gauge

Transtracheal Catheters

Transtracheal Catheters provide tracheal access for our **Manual Jet Ventilators**. Attached syringe helps to determine correct needle placement in the trachea by air aspiration.


- Packaged sterile — Ready to Use
- Unopened packages have a shelf life of a minimum of **2 years** from purchase
- Works with **BE 183-SUR & BE 183-SU Manual Jet Ventilators**
- **Foam Neck Tape with Velcro Lock**
- Attached **Syringe**
- Packaged: **Individually**



[Click to View](#)
Installation
& Usage Directions

ACU 1060.1
Transtracheal Catheter
■ 13G — Adult
[Buy Online](#)

ACU 1060.2
Transtracheal Catheter
■ 14G — Child
[Buy Online](#)



Works with **BE 183 Series** Manual Jet Ventilators **Page 20**

Oropharyngeal Airways

Oropharyngeal Airways keep the airway open by holding the tongue forward and out of the way.

- **Flange** rests against front teeth to prevent swallowing airway
- **"I Beam" Shape** allows a suction catheter to be used on either side
- Low Density Polyethylene
- Packaged: **10 or 50**
(Also available, a package of **One each of all seven SH 77 Series sizes**)



SH 77 Series
Non-Center-Vented Oropharyngeal Airways

- | | | | |
|-------------|--------------|---------|----------------------------|
| A — SH 77-1 | Neonatal | - 40mm | Buy Online |
| B — SH 77-2 | Infant | - 50mm | Buy Online |
| C — SH 77-3 | Small Child | - 60mm | Buy Online |
| D — SH 77-4 | Child | - 70mm | Buy Online |
| E — SH 77-5 | Small Adult | - 80mm | Buy Online |
| F — SH 77-6 | Medium Adult | - 90mm | Buy Online |
| G — SH 77-7 | Large Adult | - 100mm | Buy Online |



SH 77-S
One each of all seven sizes of SH 77 Series Non-Center-Vented Oropharyngeal Airways
[Buy Online](#)



We Care About the Earth

Many of our products are reusable, which helps to reduce waste! We have taken major strides to reduce manufacturing waste by recycling whenever possible. We have solar panels installed on the roof of our building to help offset our energy usage.

Laryngoscope Blades & Handles

Laryngoscope Blades

Laryngoscope Blades are available in two types, **Macintosh (curved)** and **Miller (straight)** in a variety of sizes. All blades are lighted.

- Adapts to handles with **Howland locks**
- Replacement Bulbs available
- Stainless Steel
- Packaged: **Individually**



Macintosh Blades

Curved

Laryngoscope Blades



AM MAC Series

Macintosh Blades

- A — AM MAC-1 Size 1 - 87mm
- B — AM MAC-2 Size 2 - 108mm
- C — AM MAC-3 Size 3 - 130mm
- D — AM MAC-4 Size 4 - 158mm

[Buy Online](#)

Miller Blades

Straight

Laryngoscope Blades



AM MIL Series

Miller Blades

- A — AM MIL-0 Size 0 - 75mm
- B — AM MIL-1 Size 1 - 102mm
- C — AM MIL-1.5 Size 1.5 - 127mm
- D — AM MIL-2 Size 2 - 155mm
- E — AM MIL-3 Size 3 - 195mm
- F — AM MIL-4 Size 4 - 205mm

[Buy Online](#)

Laryngoscope Handles

Laryngoscope Handles have a knurled finish for better grip.

- Fits all of our **Macintosh & Miller Blades** and adapts to most "hook-on" blades
- Stainless Steel
- Packaged: **Individually**



SM-MH
Medium Laryngoscope Handle

- Uses 2 "C" batteries

[Buy Online](#)



SM-PH
Penlite Laryngoscope Handle

- Uses 2 "AA" batteries

[Buy Online](#)



Macintosh & Miller Blade Replacement Bulbs



Replacement Bulbs for: **AM MAC-1, AM MAC-2, AM MIL-0, AM MIL-1, AM MIL-1.5, AM MIL-2, AM MIL-3**

- SM-501 — Small Frosted
- SM-503 — Small Pre-Focused

- Packaged: **Individually**

[Buy Online](#)

Replacement Bulbs for: **AM MAC-3, AM MAC-4, AM MIL-4**

- SM-502 — Large Frosted
- SM-504 — Large Pre-Focused
- SM-505 — Parabolic Reflector


- Packaged: **Individually**

[Buy Online](#)

Endotracheal Tube Changers

Endotracheal Tube Changers provide a simple technique for changing endotracheal tubes.

- No Laryngoscope needed!
- Rigid but flexible
- **Rounded tips** allow smooth, trauma-less insertion
- **Graduated marks** for easy and accurate placement of new tube
- Available in **nine** sizes for **2.5mm I.D. to 7.5mm I.D.** endotracheal tubes
- High Density Polyethylene
- Packaged: 10
(Also available, a package of **one each** of all nine JEM Series sizes)



Buy Online

Click to View

Installation & Usage Directions
Consignes d'installation et d'utilisation

JEM Series

Endotracheal Tube Changers

- A — JEM 400 fits 7.5mm I.D. tubes Length: 29"
- B — JEM 370 fits 7.0mm I.D. tubes Length: 27.5"
- C — JEM 365 fits 6.5mm I.D. tubes Length: 26.5"
- D — JEM 360 fits 6.0mm I.D. tubes Length: 26"
- E — JEM 355 fits 5.5mm I.D. tubes Length: 25"
- F — JEM 350 fits 5.0mm I.D. tubes Length: 23"
- G — JEM 340 fits 4.0mm I.D. tubes Length: 19.5"
- H — JEM 330 fits 3.0mm I.D. tubes Length: 16.5"
- I — JEM 325 fits 2.5mm I.D. tubes Length: 15"

JEM 401

One each of all nine sizes of JEM Series Endotracheal Tube Changers

Product Info

Size is clearly marked


Graduated marks for easy, accurate placement of new tube



How they work


JEM Series Endotracheal Tube Changer

1. The JEM Endotracheal Tube Changer is placed into the Endotracheal Tube. The guide mark location is noted.
2. The Endotracheal Tube is removed by passing it over the JEM Endotracheal Tube Changer. (Figure A)
3. The replacement Endotracheal Tube is inserted over the JEM Endotracheal Tube Changer. Guide marks facilitate accurate placement of the new Endotracheal Tube. (Figure B)
4. The JEM Endotracheal Tube Changer is then removed.



Pilot Tube Repair Kit

Pilot Tube Repair Kit allows cuff re-inflation on an **endotracheal** or **tracheostomy** tube with a damaged pilot line. Kit includes **Pilot Tube Repair Device** and **10 cc/mL Syringe**.



BE 409
Pilot Tube Repair Kit

- Essential for all ICU and emergency departments
- For use with endotracheal and tracheostomy tubes
- Includes **10 cc/mL Syringe** for cuff re-inflation
- Fits pilot lines of most tubes: **0.032 in. I.D. (21g)** to **0.050 in. I.D. (18g)**
- Stainless Steel Tapered Needle
- Packaged: 5 or 10







[Buy Online](#)

[Click to View](#)
Installation & Usage Directions
Consignes d'installation et d'utilisation

How it works

BE 409
Pilot Tube Repair Kit

- 1 **Pilot Tube** is cut between damaged area and patient. (Figure A)
- 2 **Pilot Tube Repair Device's Tapered Needle** is inserted into cut **Pilot Tube**. (Figure B)
- 3 **10 cc/mL Syringe** is inserted into **Pilot Tube Repair Device's Inflation Port** (Figure C), or...into **Cuff Pressure Monitor** inserted into **Pilot Tube Repair Device's Inflation Port** (Figure D).
4. Cuff is inflated gradually.



Cuff Pressure Monitors

Cuff Pressure Monitors provide a simple, effective, and economical way to avoid excessive cuff pressure.

- Packaged: **Individually**



[Click to View](#)
Installation & Usage Directions



NS 60-TBS-CP
Cuff Pressure Monitor

- Reads $\pm 60\text{cm H}_2\text{O}$ in 2cm graduations
- Bottom Mount **1/4" NPT Gauge & Adapter**
- Anodized Aluminum & Steel Adapter

[Buy Online](#)



NS 60-TRS-CP
Cuff Pressure Monitor

- Reads $\pm 60\text{cm H}_2\text{O}$ in 2cm graduations
- Rear Mount **1/4" NPT Gauge & Adapter**
- Anodized Aluminum & Steel Adapter

[Buy Online](#)



BE 148-4
Cuff Pressure Monitor Adapter

- Fits **1/4" NPT Gauges**
- Anodized Aluminum & Steel

[Buy Online](#)




NS 60-TBS-CP/NS 60-TRS-CP Gauge info for California customers:

WARNING: These products can expose you to chemicals including Lead, which is known to the State of California to cause cancer and birth defects or other reproductive harm. For more information go to www.P65Warnings.ca.gov.

How they work

NS Series Cuff Pressure Monitors (including BE 148-4)

Cuff Pressure Monitor Adapter is inserted into the **Endotracheal Tube's Pilot Line Inflation Valve**.



If cuff pressure reading is not optimal, a **10cc Syringe (not included)** is inserted into the **Cuff Pressure Monitor Adapter Valve** to increase or decrease pressure.





Swivel Barbed Tubing Adapter (Swivel "Christmas Tree" Connectors)

NEW!

Our **NEW USA-Made Swivel Barbed Tubing Adapter** helps when attaching small-bore tubing to flowmeters. The **SR 75-C** swivels **360° until tightened** while attaching. Precision-engineered and meticulously manufactured to help in delivering prescribed gas flow from flowmeters. A poorly manufactured or improperly attached **Swivel Barbed Tubing Adapter** can compromise prescribed gas flows by leaking.

- Fits standard 9/16-18 (DISS 1240) outlet
- Clear color provides universal use
- Acrylic Copolymer
- Packaged: 25 or 50



SR 75-C
Swivel Barbed Tubing Adapter

[Buy Online](#)



Barbed Tubing Adapters (“Christmas Tree” Connectors)

Barbed Tubing Adapters are precision-engineered to help in delivering prescribed gas flow from flowmeters. Poorly manufactured or improperly attached **Barbed Tubing Adapters** can compromise prescribed gas flows by leaking.

- Fits standard 9/16-18 (DISS 1240) outlet
- Four colors for easy color coding
- Packaged: 25 or 50



Make them your own!

Barbed Tubing Adapters molded with your logo, etc.

Call our friendly customer service for information.

SH 74-C
Barbed Tubing Adapter

- Clear
- Acrylic Copolymer



[Buy Online](#)

SH 74-G
Barbed Tubing Adapter

- Green
- ABS (Acrylonitrile Butadiene Styrene)



[Buy Online](#)

SH 74-W
Barbed Tubing Adapter

- White
- ABS (Acrylonitrile Butadiene Styrene)



[Buy Online](#)

SH 74-Y
Barbed Tubing Adapter

- Yellow
- ABS (Acrylonitrile Butadiene Styrene)



[Buy Online](#)

Are your current **Christmas Tree Connectors** leaking?...
test them

1. Turn flowmeter on.
2. Set flow to 3 l/min.
3. Pinch patient end of tube.

4. If ball falls to zero, it is **NOT** leaking.
 5. If ball continues to float at any level, it **IS** leaking.
- If yours are leaking, try ours!*

O₂ Instant Flow Valves

O₂ Instant Flow Valves stop the wasteful practice of oxygen dumping, while assuring the flow of oxygen to resuscitation bags. They are essential where O₂ to a **resuscitation bag** is always on, but used only when needed.

- For use with adult, pediatric, or neonatal resuscitation bags weighing more than 3 oz. (90 g)
- Flowmeter is not affected in any way!
- Standard Plated Brass DISS Fittings
- Plated Brass **Small Bore Tubing Bubble Fittings**
- Anodized Aluminum & Steel
- Packaged: **Individually**


[Click to View](#)
Installation & Usage Directions
Consignes d'installation et d'utilisation




How they work

O₂ Instant Flow Valves

When **resuscitation bag** is placed on **O₂ Instant Flow Valve** arm, oxygen flow is **stopped**.



When **resuscitation bag** is removed from **O₂ Instant Flow Valve** arm, oxygen flow **resumes**.



We did the **math...**
you do the **saving!**

Constantly flowing oxygen
without waste!



** Assuming the O₂ is not needed for various reasons 1/3 of the year. Your results may vary.*

$$\begin{aligned}
 & \$1.40/100 \text{ ft}^3 \quad \text{O}_2 \text{ cost} \\
 & \quad \times 15 \text{ lpm} \quad \text{O}_2 \text{ flow rate} \\
 & \quad \times 0.0353 \quad \text{Conversion of liters to ft.}^3 \\
 & \quad \times 175,200^* \quad \text{Minutes per year of wasted O}_2 \text{ to a resuscitation bag} \\
 & = \mathbf{\$1,299} \quad \text{Annual savings per bed!}
 \end{aligned}$$

Flowmeter-Mounted (Direct Connection to Flowmeter)



Pole-Mounted (Remote Connection to Flowmeter) *Replacement Knobs are located on pg. 12*



O₂ Instant Flow Valves *(Continued)*

Rail-Mounted (Remote Connection to Flowmeter) *Replacement Knobs are located on pg. 12*

BE 3000-RA
Rail-Mounted
O₂ Instant Flow
Valve

- Forked Arm
- Includes:
BE 1000-3
4.5-ft High Pressure
Hose Assembly
- Does not require
Bag Hanger
- Does not fit on rails flush to wall

BE 1000-3
4.5-ft. High Pressure
Hose Assembly
included with
BE 3000-RA/RS

Standard Male DISS Fitting

Small Bore Tubing Bubble Connection

Forked Arm

[Buy Online](#)

BE 3000-RS
Rail-Mounted
O₂ Instant Flow Valve

- Short Arm
- Includes:
BE 1000-3
4.5-ft High Pressure
Hose Assembly
- Includes:
BE 1000-2
Bag Hanger
- Does not fit on
rails flush to wall

Short Arm

Standard Male DISS Fitting

Small Bore Tubing Bubble Connection

BE 1000-2
Bag Hanger included
with all Short Arm
O₂ Instant Flow Valves

[Buy Online](#)

Bag Hangers

O₂ Instant Flow Valve
Bag Hangers attach resuscitation
bags to **Short Arm O₂ Instant Flow Valves**.
Not made with natural rubber latex.

BE 1000-2
Bag Hanger

- Silicone Rubber
- Packaged:
Individually

[Buy Online](#)

MADE IN USA!

R
REUSABLE

NP 1000-2D
Bag Hanger

- Neoprene
- Packaged: 50

[Buy Online](#)

MADE IN USA!

SPU
SINGLE PATIENT USE

Not made with natural rubber latex!

4.5-ft. High Pressure Hose Assembly

BE 1000-3
Replacement 4.5-ft High Pressure
Hose Assembly

- Replacement for: **BE 3000-PA & PS, BE 3000-RA & RS**
- See pg. 28 for more info
- Packaged: **Individually**

[Buy Online](#)

R
REUSABLE

Wall Mount

Wall Mount allows a Rail-Mounted O₂
Instant Flow Valve (**BE 3000-RA/RS**) to be
used in rooms without installed rails.

BE 3000-WM
Wall Mount

- Can be moved from room to room
- Anodized Aluminum
- Packaged: **Individually**

[Buy Online](#)

R
REUSABLE

MADE IN USA!

Flowmeter/Air Compressor Adapter

Flowmeter/Air Compressor Adapter
connects a flowmeter or an air compressor
directly to **22mm I.D. tubing**. Can be used
with high gas flows for CPAP.

- Standard Plated Brass Female DISS Fitting
- Security chain to prevent loss
- Anodized Aluminum & Steel
- Packaged: **Individually**



CB 23
Flowmeter/Air Compressor
Adapter

[Buy Online](#)

R
REUSABLE

MADE IN USA!

5-in-1 Connector

5-in-1 Connector helps in
connecting various sizes of
tubing or to extend the
length of tubing.

- Barbed edges help prevent disconnects
- Fits tubing with inside diameters from **3/16" to 7/16"**
- Polypropylene
- Packaged: 25 or 50



SH 7-2
5-in-1 Connector

[Buy Online](#)

SPU
SINGLE PATIENT USE

MADE IN USA!

Breathing Bags

Breathing Bags are made of neoprene, and have a 22mm I.D. connection. They can be used in closed or semi-closed (with CO₂ absorber) rebreathing circuits. **Not made with natural rubber latex.**

Three-pleat style


- Green Neoprene, Polyisoprene, & Natural Thermoplastic Rubber



KC 00.5N
.5-Liter Breathing Bag

- Packaged: **Individually**

[Buy Online](#)



SG 066-2N
1-Liter Breathing Bag

- Packaged: **1 or 25**

[Buy Online](#)



KC 002N
2-Liter Breathing Bag

- Packaged: **Individually**

[Buy Online](#)



KC 003N
3-Liter Breathing Bag

- Packaged: **Individually**

[Buy Online](#)



Two-pleat style


- Green Neoprene, Polyisoprene, & Natural Thermoplastic Rubber



KC 00.5N-2P
.5-Liter Breathing Bag

- Packaged: **Individually**


[Buy Online](#)



SG 066-2N-2P
1-Liter Breathing Bag

- Packaged: **1 or 25**

[Buy Online](#)



KC 002N-2P
2-Liter Breathing Bag

- Packaged: **Individually**

[Buy Online](#)



KC 003N-2P
3-Liter Breathing Bag

- Packaged: **Individually**

[Buy Online](#)



High Pressure Hose Assemblies

High Pressure Hose Assemblies connect oxygen remotely. They deliver a high pressure flow of oxygen from a flowmeter or any delivery source. **Not made with natural rubber latex.**

- Standard Plated Brass **Female DISS** fittings on both ends
- Medical Grade PVC Tubing
- Packaged: **Individually**



BE 1000-3
4.5-ft. High Pressure Hose Assembly

- Replacement for: **BE 3000-PA/PS & RA/RS O₂ Instant Flow Valves**

[Buy Online](#)




BE 183-6
6-ft. High Pressure Hose Assembly

- Replacement for: **BE 183 Series Manual Jet Ventilators**

[Buy Online](#)

Package Insert/User Instruction sheets for our products are available online at iiimedical.com



Flowmeter Security Devices

Flowmeter Security Devices can help eliminate the misplacement of flowmeters. They attach flowmeters securely to an oxygen wall plate.

- Packaged: **Individually**



BE 167
Flowmeter Security Wire
 ■ 8" long
 ■ Steel Components & Nylon-Coated Wire




BE 168
Flowmeter Security Chain
 ■ 8½" long
 ■ Stainless Steel Components



How they work

Flowmeter Security Devices

Attach **Flowmeter Security Device** to wall plate and Flowmeter.



O₂ Tank Wrench

O₂ Tank Wrench is made of plated steel. 12" security cable helps prevent misplacement of wrenches.

- Packaged: **Individually**



BE 169
O₂ Tank Wrench w/12" Security Cable
 ■ Steel Components & Nylon-Coated Cable



BE 169-W
O₂ Tank Wrench
 ■ Plated Steel




BE 169-C
12" Security Cable
 ■ Steel Components & Nylon-Coated Cable



How they work

O₂ Tank Wrench




Patient Accessories

Nose Clip

Nose Clip helps to hold a patient's nostrils closed while administering treatment or during clinical testing.

- Ethylene Vinyl Acetate (EVA)
- Polyoxymethylene (Acetal) (POM)
- Packaged: **50 or 100**

NC 12585
Nose Clip



Breast Marker

Breast Marker is used to assist in nipple placement during breast re-constructive surgery. A temporary impression is made by pressing the Breast Marker onto skin with a gentle rotating motion.

- Two marking edges: **38mm & 42mm**
- Methylene blue may be used as a marking agent
- Polycarbonate
- Packaged: **25**

NL 100
Breast Marker



Mask/Filter Storage Boxes

Mask/Filter Storage Boxes provide immediate availability of resuscitator bags, airways, masks, filters, eye protection, etc.

- Available in 3 sizes
- Includes 5 Red & 5 Yellow Replaceable Safety Seal Locks which ensure the integrity of box contents
- Contents clearly visible
- Can be mounted on any smooth surface
- Clear acrylic with opaque backing
- Packaged: Individually



BE 4500
Mask/Filter Storage Box

- 9¼" High
- 6" Wide
- 4½" Deep

[Buy Online](#)



BE 4500-18D
Mask/Filter Storage Box

- 18" High
- 6" Wide
- 5¼" Deep

[Buy Online](#)



BB 4500-20D
Mask/Filter Storage Box

- 20" High
- 8" Wide
- 6" Deep

[Buy Online](#)

Replacement Parts

Replacement Safety Seal Locks for all of our Mask/Filter Storage Boxes.

- Helps ensure box integrity
- For: BE 4500, BE 4500-18D, & BB 4500-20D
- Packaged: 5



BE 4501
Replacement Safety Seal Locks

- Red

[Buy Online](#)

BE 4502
Replacement Safety Seal Locks

- Yellow

[Buy Online](#)



Anesthesia Flush Plug

Anesthesia Flush Plug provides a simple way to plug an anesthesia machine breathing circuit for flushing or other maintenance.

AM 100P
Anesthesia Flush Plug

- Can help reduce the possibility of exposure of personnel to anesthetic gases
- Protects anesthesia machine or mask from dust and debris when not in use
- Keeps ventilator circuits clean when not in use
- High Density Polyethylene
- Packaged: 25 or 50



[Buy Online](#)



15mm O.D.


[Click to View](#)
Installation
& Usage Directions

How it works

AM 100P
Anesthesia Flush Plug

Typical Anesthesia Machine Breathing Circuit with Gas Sampling Elbow and Anesthesia Mask.

Anesthesia Flush Plug can be used in multiple circuit locations.



"NO SMOKING - OXYGEN IN USE" Sign

"No Smoking - Oxygen in Use" Sign is sold in multiples for use throughout your facility.

BE 352

"NO SMOKING - OXYGEN IN USE" Sign

- 5" x 8"
- Red on white for high visibility
- Inexpensive way to display "No Smoking" sign where oxygen is in use
- English one side / Spanish one side
- White Vinyl
- Packaged: 50 or 100



Magnetic "OXYGEN - NO SMOKING" Signs

Magnetic "Oxygen - No Smoking" Signs are sold in multiples for use throughout your facility.

- 1¾" x 8"
- Red on white for high visibility
- Available in English or Spanish
- Packaged: 10



BE 354
Magnetic
"OXYGEN —
NO SMOKING"
Sign

- English



BE 355
Magnetic
"OXIGENO —
NO FUMAR"
Sign

- Spanish



SeCUREacord Cord Covers

SeCUREacord Cord Covers are a necessity within a healthcare setting. They are intended to help secure the **monitoring cords** attached to a patient during **movement** or **transfer**. The use of **SeCUREacord Cord Covers** helps prevent medical cords from tangling, which can cause accidental disconnection from **patient monitoring devices**. The **Cord Covers** also reduce the risk to **staff, visitors, and patients** due to trip hazards associated with uncontrolled cords. Additionally, the **SeCUREacord Cord Covers** help prevent the cords from being contaminated with bodily fluids.

- Clear tube with black Velcro enclosures on both ends
- Corrugated spiral design to prevent kinking
- Approximate 1" diameter
- Available in 6 sizes
- Made of Polyolefin Plastomer & Nylon
- Sold: **Individually**

SJ SERIES

SeCUREacord Cord Covers

- | | |
|---------------------|---------------------|
| SJ 118 — 18" length | SJ 136 — 36" length |
| SJ 124 — 24" length | SJ 142 — 42" length |
| SJ 130 — 30" length | SJ 148 — 48" length |



NEW!

[Click to View Installation & Usage Directions](#)

MADE IN USA!

SPU SINGLE PATIENT USE



Product Index

A	BE 142.....10	BE 4501.....30	N	SG 066-2N.....28
ACU 1060.1.....21	BE 145-1522P.....7	BE 4502.....30	NB 105-6.....5	SG 066-2N-2P.....28
ACU 1060.2.....21	BE 145-2230P.....7	C	NC 12585.....29	SG 101.....5
ADP 150.....7	BE 145-2P.....7	CAP-S.....8	NL 100.....29	SG 101-D.....5
AM 100P.....30	BE 147-GC.....19	CB 23.....27	NP 1000-2D.....27	SG 1010.....4
AMMAC Series.....22	BE 147-GCV.....19	D	NP 105-1.....13	SG 102.....5
AMMIL Series.....22	BE 148-1.....18	DS 7-45.....9	NP 130-22.....6	SH 7-2.....27
B	BE 148-2.....19	E	NS 120-TRR.....17	SH 74 Series.....25
BB 4500-20D.....30	BE 148-4.....24	EH 2.....17	NS 120-TRS.....16	SH 77 Series.....21
BE 1000-2.....27	BE 148-5-1.....18	H	NS 120P-TRS.....16	SJ 100 Series.....31
BE 1000-3.....27, 28	BE 148-5-2.....18	HG 100.....15	NS 30-115-5.....10	SM-501.....22
BE 101-00-1426.....7	BE 148-5-2.5.....18	J	NS 30-115-5L.....10	SM-502.....22
BE 101-00-3443.....7	BE 148-6.....18	JEM Series.....23	NS 30-115-U.....10	SM-503.....22
BE 102.....13	BE 149-1.....18	K	NS 30-115-UL.....10	SM-504.....22
BE 105-101-40.....6	BE 150-5G.....11	KC 00.5N.....28	NS 30-PBR.....17	SM-505.....22
BE 105-2222.....3	BE 150-5LF.....11	KC 00.5N-2P.....28	NS 30-PBS.....16	SM-MH.....22
BE 105-2X22.....4	BE 16-2.....4	KC 002N.....28	NS 30-TRR.....17	SM-PH.....22
BE 105-2X22D.....4	BE 167.....29	KC 002N-2P.....28	NS 30-TRS.....16	SR 75-C.....25
BE 105-40.....6	BE 168.....29	KC 003N.....28	NS 60-PBR.....17	W
BE 105-7.....9	BE 169.....29	KC 003N-2P.....28	NS 60-PBS.....16	WT 125.....11
BE 105-7-1.....9	BE 169-C.....29	KC 003N-2P.....28	NS 60-TBR.....17	WT 125-BL.....11
BE 105-8.....9	BE 169-W.....29	KC 1.....12	NS 60-TBS.....16	WT 125-BOT.....11
BE 106.....3	BE 183-2.....20	KC 10.....13	NS 60-TBS-CP.....24	WT 250.....11
BE 114-5.....12	BE 183-6.....20, 28	KC 1010.....4	NS 60-TRR.....17	WT 250-BL.....11
BE 114-K.....12	BE 183-G.....20	KC 105.....4	NS 60-TRS.....16	WT 250-BOT.....11
BE 114-5K.....12	BE 183-R.....20	KC 145-1522D.....7	NS 60-TRS-CP.....24	
BE 114-5TB.....13	BE 183-SU.....20	KC 145-2230D.....7	NS-GC.....19	
BE 115-2N.....10	BE 183-SUR.....20	KC 1515.....4	NS-GS.....19	
BE 117.....9	BE 190.....5	KC 18901-K.....6	NS-RC.....17	
BE 117-W.....4	BE 191.....5	KC 18901-L.....6	NS-SC.....16	
BE 120-22D.....3	BE 192.....5	KC 18901-O.....6	P	
BE 120-AKW.....3	BE 193.....5	KC 18901-P.....6	PF 001.....7	
BE 122-24.....12	BE 195-1522A.....5	KC 211-D.....8	PF 1010.....7	
BE 122-2TB.....13	BE 215-2A.....10	KC 2119.....5	PF 1510.....7	
BE 122-2TSB.....13	BE 30-115-B.....10	KC 212-D.....8	PF 1515.....7	
BE 122-30.....12	BE 30-115-BL.....10	KC 2124.....4	R	
BE 122-S.....12	BE 30-142.....10	KC 2126.....5	RB 002.....12	
BE 125.....4	BE 3000-FA.....26	KC 213-D.....8	RB KNOB.....12	
BE 125-1.....4	BE 3000-FS.....26	KC 5.....6	RTC 15-D.....3	
BE 125-1T.....4	BE 3000-PA.....26	KC 6PA.....6	RTC 15-P.....3	
BE 125-1TD.....4	BE 3000-PS.....26	KC Series.....14	RTC 17-1.....3,9	
BE 127.....8	BE 3000-RA.....27	KC-N Series.....14	RTC 21-A.....1	
BE 127-40.....8	BE 3000-RS.....27	KC-P Series.....14	RTC 21-D.....1	
BE 127-7.5.....8	BE 3000-WM.....27	KC-3 Series.....14	RTC 22-A.....3	
BE 127-7.5-40.....8	BE 352.....31	KC-5 Series.....15	RTC 22-D.....3	
BE 130.....11	BE 354.....31	KC-5N Series.....15	RTC 23-A.....3	
BE 130-22B.....11	BE 355.....31	KC-5P Series.....15	RTC 23-D.....3	
BE 130-23BB.....11	BE 409.....24	KM 030.....9	RTC 24-V.....2	
BE 131.....11	BE 40A.....6	M	RTC 24-V KIT.....2	
BE 139.....6	BE 40H.....6	MS 101-40.....6	RTC 24-VP.....2	
BE 139-105.....5	BE 4500.....30		RTC 26-C.....1	
	BE 4500-18D.....30			

