

You can now use
**New Generation
 MDI Canisters** with
 a dose counter inline
 on a patient circuit!




Introducing...

Instrumentation Industries, Inc.
 newest MDI Adapter — **RTC 24-V**

You can now easily use the **New Generation MDI Canisters** with a dose counter on ventilator patients! No more prying, breaking, etc., of dose counters!


The **RTC 24-V** also accepts **Standard MDI Canisters** for added versatility.

Now you can administer **Bronchodilators, Corticosteroids, or other Aerosolized Medications** from either type of MDI Canister inline on patient circuits.



RTC 24-V

- Inlet 22mm I.D.
- Outlet 15mm I.D./22mm O.D.
- Clearly marked Medication Flow Direction Arrows
- Fits either New Generation MDI Canisters with dose counters or Standard MDI Canisters
- Counter Activator
- Thermoplastic Rubber Cap
- Sold in a package of 25 or 50
- Single Patient Reuse*



No need to pry or break off dose counters on New Generation MDI Canisters!



Also available... RTC 24-V KIT

RTC 24-V KIT allows for versatility in creating a reservoir for the **RTC 24-V MDI Adapter** either at the endotracheal tube, **1** or in the inspiratory limb. **2**

RTC 24-V KIT includes a **RTC 24-V MDI Adapter** along with our **KC 2124 15mm I.D./22 O.D. x 15mm I.D./22mm O.D. Adapter** and our **KC 004-5 4" length of Clear Tubing**.



- **RTC 24-V**
- Inlet 22mm I.D.
- Outlet 15mm I.D./22mm O.D.
- **KC 2124**
- 15mm I.D./22mm O.D. x 15mm I.D./22mm O.D.
- **KC 004-5**
- 22mm I.D. x 22mm I.D.
- 4" Length
- **RTC 24-V KIT**
- Sold in a package of 10
- Single Patient Reuse*



ways to order!

Order online
www.iiimedical.com

Order by phone
 1-800-633-8577

*Single Patient Reuse — Can be used to repeatedly administer medications to one single patient while installed in a breathing circuit.




Instrumentation Industries, Inc.



Reliable patient circuit products!

Patient Wye Connectors

Reusable or Single Patient Use


 **Patient Wye Connectors** have industry standard dimensions for adult circuits. They're available in various types — reusable, single patient use, 60°, 90°, with or without pressure, temperature or nebulizer ports.

Visit www.iiimedical.com or all of our Patient Circuit Products.

Single Patient Use →
sold in packages of 50 & 100

Reusable →
sold in packages of 1 & 10

Single Patient Use
22mm O.D. x 22mm O.D. x 22mm O.D.



Reusable — 22mm O.D. x 22mm O.D. x 15mm I.D./22mm O.D.



Ventilator Tubing

Mechanical ventilation is a life-saving procedure.

 **Ventilator Tubing** is corrugated to prevent kinking. It has a smooth internal diameter which promotes laminar gas flow, and resists the trapping of condensation. All of our ventilator tubing has industry standard connection diameters of 10, 15, & 22mm I.D.

Stock lengths are sold individually...No minimum!


Using securely attached tubing is extremely important.



Custom lengths are available (minimum orders are required) call for more info

2 ways to order!

 Order online
www.iiimedical.com

 Order by phone
1-800-633-8577

 **Instrumentation Industries, Inc.**



 **Instrumentation Industries, Inc.**

2990 Industrial Boulevard • Bethel Park, PA 15102

Address Service Requested

Mailroom Please Route to:

Respiratory Care/Cardiopulmonary

or Anesthesia Critical Care Emergency Care

PRST STD
US POSTAGE
PAID
PITTSBURGH, PA
PERMIT NO 1529

Plastic
Metals
Silicone Rubber

 **MedicalOEM**
A division of Instrumentation Industries, Inc.

Do you have an idea for a Respiratory Care Product?

Our **Medical OEM Division** can develop and manufacture your products. From single prototypes to full production runs, we have the facility and staff to meet your needs!



Instrumentation Industries, Inc.
is an ISO 13485: 2003 certified company.

iii1209

Finally, the problem is solved!

look inside...

 **Instrumentation Industries, Inc.**

 **Order Online!!**
www.iiimedical.com

