











**Instrumentation  
Industries, Inc.**  
*Since 1967*

# Packaging Symbols

The following is a glossary of symbols used on Instrumentation Industries, Inc. product labeling.

SYMBOLS GLOSSARY			
The following Symbols have been adopted into, and are listed in International Standard ISO 15223-1:2012.			
Reference Symbol	Title of Symbol	Description of Symbol	ISO 7000 Reg. no.
	Manufacturer	Indicates the medical device manufacturer, as defined in EU Directives 90/385/EEC, 93/42/EEC and 98/79/EC.	3082
	Non-sterile	Indicates a medical device that has not been subjected to a sterilization process.	2609
	Do not re-use	Indicates a medical device that is intended for one use, or for use on a single patient during a single procedure.	1051
	Consult instructions for use	Indicates the need for the user to consult the instructions for use.	1641
	Catalogue number	Indicates the manufacturer's catalogue number so that a specific medical device can be identified.	2493
	Batch code	Indicates the manufacturer's batch code so that the batch or lot can be identified.	2492
The following Symbols are not listed in International Standard ISO 15223-1:2012, but are included on product labeling, when applicable.			
	Rx Only	Federal law restricts this device to sale by or on the order of a physician.	N/A
	MR Unsafe	Indicates a medical device poses unacceptable risks to the patient, medical staff or other persons within the MR environment.	N/A