

Why use our Gauges?

We have 6 reasons...



#1 Accuracy

Vacuum/Pressure Gauges


Vacuum/Pressure Gauges provide accurate pressure monitoring in incremental measurements as low as 1cm H₂O. Needle can easily be "re-zeroed" for accuracy.

- ◆ Reusable
- ◆ Sold Individually

Vacuum/Pressure Gauges		
a	NS 30-PBS NS 60-PBS	±30cm H ₂ O in 1cm increments Bottom 1/8" Small Bore Tubing ±60cm H ₂ O in 2cm increments Bottom 1/8" Small Bore Tubing
b	NS 30-TRS NS 60-TRS NS 120-TRS	±30cm H ₂ O in 1cm increments Rear 1/4" NPT ±60cm H ₂ O in 2cm increments Rear 1/4" NPT ±120cm H ₂ O in 4cm increments Rear 1/4" NPT
c	NS 60-TBS	±60cm H ₂ O in 2cm increments Bottom 1/4" NPT

Replacement crystals are available for all Vacuum/Pressure Gauges.

Scan for more info on our line of Vacuum/Pressure Gauges



#2 Quality

NIF Meters

NIF Meters capture a patient's best effort with a resettable pointer. Available in incremental measurements as low as 1cm H₂O. Needle can easily be "re-zeroed" for accuracy.


- ◆ Reusable
- ◆ Sold Individually

NIF Meters		
a	NS 30-PBR NS 60-PBR	±30cm H ₂ O in 1cm increments Bottom 1/8" Small Bore Tubing ±60cm H ₂ O in 2cm increments Bottom 1/8" Small Bore Tubing
b	NS 30-TRR NS 60-TRR NS 120-TRR	±30cm H ₂ O in 1cm increments Rear 1/4" NPT ±60cm H ₂ O in 2cm increments Rear 1/4" NPT ±120cm H ₂ O in 4cm increments Rear 1/4" NPT
c	NS 60-TBR	±60cm H ₂ O in 2cm increments Bottom 1/4" NPT

1/8" Small Bore Tubing connection models include BE 148-6 Gauge Adapter.

1/4" NPT connection models include BE 149-1 Gauge Adapter. Replacement crystals are available for all NIF Meters.

Scan for more info on our line of NIF Meters



Small Bore Medical Tubing

Small Bore Medical Tubing meets USP Class VI, ISO 10993, and REACH requirements. Made of Tygon® ND 100-65.

1/8" I.D.		3-16" I.D.	
SB 1/8-01	1-ft.	SB 3/16-01	1-ft.
SB 1/8-06	6-ft.	SB 3/16-06	6-ft.
SB 1/8-25	25-ft.	SB 3/16-25	25-ft.



Pressure Adapters

Pressure Adapters can be placed in any adult, pediatric, or neonatal ventilator circuit. Attach one of our Vacuum/Pressure Gauges or NIF Meters for inline pressure measuring. 7 models available with various port styles. Visit iiimedical.com for more info.



#3 Convenience

Gauge Clamps

Gauge Clamps offer the convenience of mounting Vacuum Pressure Gauges and NIF Meters in a variety of places. Mount gauges with threaded 1/4" NPT connections to poles, rails, chairs, stretchers, incubators, ventilators, etc. Precision machined by skilled craftsmen from sturdy aluminum and then anodized.

- ◆ Reusable
- ◆ Sold Individually

Gauge Clamps		
Gauge Pole Clamps		
a	BE 148-5-1	Fits 3/4" diameter poles
b	BE 148-5-2	Fits 3/4" to 1 1/2" diameter poles
Gauge Clamps		
c	BE 148-5-2.5	Clamp opens 1" to 2 1/2"
d	BE 148-5-3.5	Clamp opens 2 1/2" to 3 1/2"
Rail Clamp		
e	RB 004	Fits Fairfield-Style Rails



Scan for more info on our line of Gauge Clamps

#4 Adaptability

Gauge Adapters

Gauge Adapters allow for easy remote inspiratory measurement...or change 1/4" NPT to 1/8" small bore tubing connections. Designed for all Instrumentation Industries, Inc. Vacuum/Pressure Gauges and NIF Meters.

- ◆ Reusable
- ◆ Sold Individually

Gauge Adapters		Vent	
a	BE 148-1	1/4" NPT x 1/8" small bore	No
b	BE 148-6	15mm I.D./22mm O.D. x 1/8" small bore	Yes
c	BE 149-1	15mm I.D./22mm O.D. x 1/4" NPT	Yes



Scan for more info on our line of Gauge Adapters

#5 Assurance

Gauge Cover

Gauge Cover assures you that our Vacuum/Pressure Gauges and NIF Meters will have the highest level of protection from falls, drops, vibrations, and other shocks. Unique "clam shell" design securely locks the gauge in place. The ribbed design provides maximum shock protection. Made of thermoplastic rubber*. Designed to fit all of our NS Series NIF Meters and Vacuum/Pressure Gauges.

- ◆ Reusable
- ◆ Sold Individually

Gauge Cover	
NS-GC	

*Not manufactured with natural rubber latex.



Scan for more info on our line of Gauge Covers

#6 Stability

Gauge Stand

Gauge Stand keeps a gauge upright and easy to view. "No-skid" rubber feet keep the whole unit stable. Made of lightweight polystyrene-butadiene, it can be placed on a variety of surfaces. Specially designed for our NS Series -PBR, -PBS, -TRR, & -TRS Vacuum/Pressure Gauges and NIF Meters.

- ◆ Reusable
- ◆ Sold Individually

Gauge Stand	
NS-GS	

Gauges are easily & quickly inserted!



Scan for more info on our Gauge Stand

3 ways to order!

online iiimedical.com

phone
1-800-633-8577 (inside USA)
1-412-854-1133 (outside USA)

fax
1-877-633-8661 (inside USA)
1-412-854-5668 (outside USA)

We care about the



Earth!
Many of our products are reusable which can help to reduce waste

Instrumentation Industries, Inc.
Since 1967






There when you need it!

Facility Accessories

Mask/Filter Storage Boxes

Mask/Filter Storage Boxes provide immediate availability of resuscitator bags, airways, masks, filters, eye protection, etc. Made of clear acrylic for easy visibility. Can be mounted on any smooth surface. Now available in 3 sizes. All three sizes include 5 Red & 5 Yellow Replaceable Seal Locks to ensure the integrity of box contents. Replacement Seal Locks are available.

- ◆ Reusable
- ◆ Sold Individually



Mask/Filter Storage Boxes				
a	BE 4500	Height 9 1/4"	Width 6"	Depth 4 1/2"
b	BE 4500-18D	Height 18"	Width 6"	Depth 5 1/4"
c	BB 4500-20D	Height 20"	Width 8"	Depth 6"

Scan for more info on our line of Mask/Filter Storage Boxes




Replacement Seal Locks

Replacement Seal Locks for all of our Mask/Filter Storage Boxes. Available in Red & Yellow.

- ◆ Single Use
- ◆ Sold in Packages of 5

Seal Locks		
a	BE 4501	Red
b	BE 4502	Yellow



Oxygen/No Smoking Signs


Oxygen/No Smoking Signs are essential where oxygen is in use. 5" high x 8" wide. Printed red on white for high visibility. Warning is printed in English on one side and Spanish on the other. Sold in multiples for use throughout your facility.

- ◆ Reusable
- ◆ Sold in Packages of 50 & 100

Oxygen/No Smoking Signs
BE 352



Scan for more info on our Oxygen/No Smoking Signs




Magnetic Oxygen/No Smoking Signs

Magnetic Oxygen/No Smoking Signs are essential where oxygen is in use. 1 3/4" wide x 8" high. Printed red on white for high visibility. Sold in multiples for use throughout your facility. Available in English and Spanish.

- ◆ Reusable
- ◆ Sold in Packages of 10

Magnetic Oxygen/No Smoking Signs				
a	BE 354	English	Height 8"	Width 1 3/4"
b	BE 355	Spanish	Height 8"	Width 1 3/4"



Breathing Bags

Breathing Bags can greatly aid in providing increased O2 concentrations. Available in 4 sizes. — .5, 1, 2, & 3 liter. Visit iiimedical.com for more info.



PRST STD
US POSTAGE
PAID
PITTSBURGH, PA
PERMIT NO 1529

Instrumentation Industries, Inc.
2990 Industrial Boulevard • Bethel Park, PA 15102

Address Service Requested

Mailroom Please Route to:

Respiratory Care/Cardiopulmonary or Anesthesia — Critical Care — Emergency Care

Do you have an idea for a Respiratory Care Product?

Plastic

Silicone Rubber

Metal

Our Medical OEM Division can develop and manufacture your products. From single prototypes to full production runs, we have the facility and staff to meet your needs!



A division of Instrumentation Industries, Inc.



Instrumentation Industries, Inc. is an ISO 13485:2003 certified company.

iii0815



Order Online!!
www.iiimedical.com

Featured Products

Vacuum/Pressure Gauges
NIF Meters
Gauge Accessories
Facility Accessories



Instrumentation Industries, Inc.
Since 1967

Accurate Pressure Monitoring

