


We're your connection!

When you need a connector or adapter for a patient circuit...we're your connection.



Adapters vs Connectors	
Adapters generally connect items of different sizes.	Connectors generally connect items of the same size.

General Adapters

General Adapters are the same size on both ends for circuit customization when needed. Available in adult, pediatric, and neonatal sizes.

Made in USA!



KC 105

General Adapters		
KC 105	Adult R	15mm I.D./22mm O.D. x 15mm I.D./22mm O.D. Packaged: 1 or 10
KC 2124	Adult SPU	15mm I.D./22mm O.D. x 15mm I.D./22mm O.D. Packaged: 25 or 100
BE 105-2X22	Adult R	22mm I.D. x 22mm I.D. Packaged: 5 or 25
BE 105-2X22D	Adult SPU	22mm I.D. x 22mm I.D. Packaged: 20 or 100
KC 1515	Pediatric SPU	15mm O.D. x 15mm O.D. Packaged: 25 or 100
KC 1010	Neonatal SPU	10mm O.D. x 10mm O.D. Packaged: 25 or 100
SG 1010	Neonatal R	10mm O.D. x 10mm O.D. Packaged: 1 or 10

Tee Adapters

Tee Adapters help add versatility to patient circuits. Available in seven types to allow for various configurations during patient circuit assembly.

Made in USA!



BE 125

Tee Adapters		
BE 125	R Packaged: 1 or 10	15mm O.D./22mm O.D. x 22mm O.D. x 15mm I.D./22mm O.D.
BE 125-1	R Packaged: Individually	22mm O.D. x 22mm O.D. x 15mm I.D./22mm O.D.
BE 125-1T	R Packaged: 1 or 10	22mm O.D. x 22mm I.D. x 15mm I.D./22mm O.D.
BE 125-1TD	SPU Packaged: 10, 20 or 100	22mm I.D. x 22mm O.D. x 15mm I.D.
BE 117-W	SPU Packaged: 10	22mm I.D. x 22mm O.D. x 15mm I.D./22mm O.D.
BE 16-2	R Packaged: Individually	22mm I.D. x 18mm I.D./22mm O.D. x 18mm I.D./22mm O.D.
SG 218	R Packaged: 1 or 10	18mm I.D./22mm O.D. x 22mm O.D. x 18mm I.D./22mm O.D.


Inline MDI Adapters

Inline MDI Adapters direct aerosol medication mainstream from Standard MDI Canisters for ventilator patients. Clearly marked with airflow direction arrows. Plastic or Metal-tipped Standard MDI Canisters may be used.

Made in USA!



RTC 22-A




Inline MDI Adapters		
RTC 15-D	SPU Packaged: 50 or 100	Inlet: 15mm O.D. Outlet: 15mm I.D.
RTC 15-P	SPU Packaged: 50 or 100	Inlet: 15mm I.D. Outlet: 15mm I.D.
RTC 22-A	R Packaged: 5 or 10	Inlet: 22mm I.D. Outlet: 15mm I.D./22mm O.D.
RTC 22-D	SPU Packaged: 50 or 100	Inlet: 22mm I.D. Outlet: 15mm I.D./22mm O.D.
RTC 23-A	R Packaged: 5 or 10	Inlet: 15mm I.D./22mm O.D. Outlet: 15mm I.D./22mm O.D.
RTC 23-D	SPU Packaged: 50 or 100	Inlet: 15mm I.D./22mm O.D. Outlet: 15mm I.D./22mm O.D.

Inline Aerosol Adapter

Inline Aerosol Adapter allows a Combivent® Respimat® Inhaler to be adapted to a ventilator circuit. Immediately before the ventilator cycle, activate the Combivent® Respimat® Inhaler. The medication then becomes suspended in the Inline Aerosol Adapter until the inspiratory cycle of the ventilator carries it to the patient. Patent Pending.

Made in USA!



Inline Aerosol Adapter		
RTC 26-C	SPR** Packaged: 25 or 50	Inlet: 22mm O.D. Outlet: 22mm I.D.

Combivent® is a registered trademark of Boehringer Ingelheim Pharmaceuticals, Inc. Copyright© 2016 Boehringer Ingelheim International GmbH. Respimat® is a registered trademark and used under license from Boehringer Ingelheim International GmbH.

Specialty Adapters

Specialty Adapters increase treatment delivery options. Choose from Adapter with a vent, Swivel Adapter for CPAP/BiPAP systems, 45° Elbow Adapter, & Mouthpiece Adapter.

Made in USA!



BE 106

Specialty Adapters		
BE 106	Adapter w/Vent R	15mm I.D./22mm O.D. x 15mm I.D./22mm O.D. Packaged: Individually
BE 105-2222	Swivel Adapter SPU	22mm I.D. x 15mm I.D./22mm O.D. Packaged: 10 or 100
BE 120-22D	Mouthpiece SPU	22mm I.D. Packaged: 100
BE 120-AKW	45° Angle Adapter SPU	22mm O.D. x 22mm I.D. Packaged: 100

Reducing Adapters

Reducing Adapters are available in 13 types! Reduce 30mm to 22mm adult size fittings, 22mm adult size fittings to 15mm pediatric size fittings, etc.

Made in USA!



Reducing Adapters			
KC 2119	SPU	Pkg: 25	15mm I.D./22mm O.D. x 15mm O.D.
KC 2126	SPU	Pkg: 25 or 100	22mm I.D. x 22mm O.D.
NB 105-6	R	Pkg: 1	15mm I.D./22mm O.D. x 22mm I.D.
BE 139-105	R	Pkg: 1 or 10	15mm I.D./22mm O.D. x 15mm O.D.
BE 190	R	Pkg: 1 or 10	30mm O.D. x 22mm O.D.
BE 191	R	Pkg: 1 or 10	30mm I.D. x 22mm I.D.
BE 192	R	Pkg: 1 or 10	30mm O.D. x 22mm I.D.
BE 193	R	Pkg: 1 or 10	30mm I.D. x 22mm O.D.
BE 195-1522A	R	Pkg: 1 or 10	30mm I.D. x 15mm I.D./22mm O.D.
SG 193-M	R	Pkg: 1 or 10	30mm I.D. x 22mm O.D.
SG 101	R	Pkg: 1 or 10	22mm I.D. x 15mm O.D.
SG 101-D	SPU	Pkg: 25 or 100	22mm I.D. x 15mm O.D.
SG 102	R	Pkg: 1 or 10	22mm I.D. x 10mm O.D.

Pressure Adapters

Pressure Adapters can be placed in any adult, pediatric, or neonatal ventilator circuit. Attach one of our Vacuum/Pressure Gauges for inline pressure measuring.

Made in USA!



BE 139

Pressure Adapters		
BE 139	Barbed Port R	15mm I.D./22mm O.D. x 15mm O.D. Packaged: 1 or 10
BE 105-40	Capped Barbed Port R	15mm I.D./22mm O.D. x 15mm I.D./22mm O.D. Packaged: Individually
BE 105-101-40	Capped Barbed Port R	15mm I.D./22mm O.D. x 22mm I.D. Packaged: Individually
KC 5	Luer Lock Port R	15mm I.D./22mm O.D. x 15mm O.D. Packaged: Individually
KC 6PA	Capped Luer Lock Port SPU	15mm I.D. x 15mm I.D. Packaged: 25 or 100
NP 130-22	Capped Barbed Port R	15mm I.D./22mm O.D. x 22mm I.D. Packaged: Individually
MS 101-40	Barbed Port R	15mm I.D. x 15mm I.D. Packaged: Individually

Universal Inline MDI Adapters


Universal Inline MDI Adapters direct aerosol medication mainstream from a Canister with a Dose Counter or a Standard MDI Canister for ventilator patients. Clearly marked with airflow direction arrows. Plastic or metal-tipped Standard Medication Canisters may be used.

Made in USA!



RTC 24-V

Universal Inline MDI Adapters		
RTC 24-V	SPR* Packaged: 25 or 50	Inlet: 22mm I.D. Outlet: 15mm I.D./22mm O.D.
RTC 24-VP	SPR* Packaged: 25 or 50	Inlet: 15mm O.D. Outlet: 15mm I.D./22mm O.D.



R	Reusable	SPU	Single Patient Use
SPR*	Single Patient Reuse (Can be used to deliver medication multiple times to a single patient while installed in a breathing circuit.)		
SPR**	Single Patient Reuse (Intended to be used to deliver Combivent® (ipratropium bromide and albuterol) Respimat® medication multiple times to a single patient while installed in a breathing circuit.)		

3 ways to order!

Online
iimedical.com

Phone 1-800-633-8577
USA toll free
1-412-854-1133
Outside USA

Fax 1-877-633-8661
USA toll free
1-412-854-5668
Outside USA



Instrumentation Industries, Inc.
Since 1967





We're your connection!

When you need a connector or adapter for a patient circuit
...we're your connection.

Patient Wye Connectors


Reusable Patient Wye Connectors are available with pressure ports and/or temperature probe ports for easy monitoring of the circuit at the patient wye.



Reusable Patient Wye Connectors	
15mm I.D./22mm O.D. x 22mm O.D. x 22mm O.D. — Packaged: 1 or 10	
BE 127	No Monitoring Ports
BE 127-40	Barbed Pressure Port
BE 127-7.5	7.5mm I.D. Temperature Monitoring Port
BE 127-7.5-40	7.5mm I.D. Temperature Probe Port & Barbed Pressure Port

Single Patient Use Patient Wye Connectors

are available with ports at the patient wye.



Single Patient Use Patient Wye Connectors		
22mm O.D. x 22mm O.D. x 22mm O.D. — Packaged: 50 or 100		
KC 211-D	60° Wye	18mm I.D./22mm O.D. Port
KC 212-D	90° Wye	18mm I.D./22mm O.D. Port
KC 213-D	60° Wye	No Additional Ports

Silicone Rubber Connectors

Silicone Rubber Connectors are just the thing to connect circuit fittings when they don't quite fit. Dissimilar sizes, or accessories of the same size, can be used together. Make use of "odd" size accessories. Autoclavable Silicone Rubber Connectors remain pliable following sterilization.



Made in USA!

Silicone Rubber Connectors		
BE 101-00-1426	Adult Tapered R	27mm I.D. x 22mm I.D. large end also fits flanged 22mm O.D. Packaged: 10 or 50
BE 101-00-3443	Adult Straight R	22mm I.D. x 22mm I.D. Packaged: 10 or 50
PF 001	Adult to Pediatric Reducing R	22mm I.D. x 15mm I.D. Packaged: 10 or 50
PF 1510	Pediatric to Neonatal Reducing R	15mm I.D. x 10mm I.D. Packaged: 10 or 50
PF 1515	Pediatric Straight R	15mm I.D. x 15mm I.D. Packaged: 10 or 50
PF 1010	Neonatal Straight R	10mm I.D. x 10mm I.D. Packaged: 10 or 50

90° Elbow Connectors

90° Elbow Connectors are available with a luer lock port or a barbed port to allow for gas sampling or pressure monitoring.



Made in USA!

90° Elbow Connectors		
KC 18901-K	No Monitoring Ports SPU	15mm I.D./22mm O.D. x 15mm O.D. Packaged: 50
KC 18901-L	Capped Luer Lock Port R	15mm I.D./22mm O.D. x 15mm O.D. Packaged: 10 or 20
KC 18901-O	Barbed Port R	15mm I.D./22mm O.D. x 15mm O.D. Packaged: 10 or 20
KC 18901-P	No Ports R	15mm I.D./22mm O.D. x 15mm O.D. Packaged: 10 or 20

Hose Connectors

Hose Connectors change a 15mm, 22mm or 30mm I.D. low pressure gas outlet source to a 1/4" small bore tubing connection. Also available, a 15mm O.D. x 1/2" O.D. hose barb connector that fits 10mm I.D. tubing.



Hose Connectors			
KC 145-1522D			
ADP 150	R	15mm O.D. x 1/2" O.D.	Pkg: 10 or 50
BE 145-2P	R	15mm O.D. x Tapered Barbed Port	Pkg: 10 or 50
KC 145-1522D	SPU	15mm I.D./22mm O.D. x Tapered Barbed Port	Pkg: 25 or 100
BE 145-1522P	R	15mm I.D./22mm O.D. x Tapered Barbed Port	Pkg: 10 or 50
BE 145-2230P	R	22mm I.D./30mm O.D. x Tapered Barbed Port	Pkg: 10 or 50

3 ways to order!



Online

iiimedical.com

Phone

1-800-633-8577
USA toll free
1-412-854-1133
Outside USA

Fax

1-877-633-8661
USA toll free
1-412-854-5668
Outside USA



Instrumentation Industries, Inc.
Since 1967



Instrumentation Industries, Inc.

2990 Industrial Boulevard • Bethel Park, PA 15102

Address Service Requested

Mailroom Please Route to:

Respiratory Care/Cardiopulmonary
or Anesthesia — Critical Care — Emergency Care

PRST STD
US POSTAGE
PAID
PITTSBURGH, PA
PERMIT NO 1529

Do you have an idea for a Respiratory Care Product?

Plastic

Silicone Rubber

Metal

Our **Medical OEM Division** can develop and manufacture your products. From single prototypes to full production runs, we have the facility and staff to meet your needs!



A division of Instrumentation Industries, Inc.



Instrumentation Industries, Inc. is an ISO 13485:2003 certified company.

iii0318

We're your connection!

When you need a connector or adapter for a patient circuit
...we're your connection.



Instrumentation Industries, Inc.
Since 1967



Order Online
iiimedical.com

